ASSKEA
_____ medical

DE

EN

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www.asskea.de

Copyright © 2016 ASSKEA GmbH, Gebesee.

The safety of the ASSKEA secretion aspirator proVisio® M28 complies with the acknowledged rules of technology and meets the requirements of the **German Medical Devices Act**.

The ASSKEA secretion aspirator proVisio® M28 bears the **CE marking CE1434** in accordance with the EU Council Directive 93/42/EEC concerning medical devices and meets the essential requirements of Annex I of this directive.

The ASSKEA secretion aspirator proVisio® M28 has been tested in accordance with IEC 62353.

The **quality management system** applied by ASSKEA GmbH is certified in compliance with the relevant international standards.

The ASSKEA secretion aspirator proVisio® M28 is a medical aspirator classified as class IIa in accordance with EU Council Directive 93/42/EEC, Annex IX.

Errors and omissions excepted.

Contents

EN

44

1	User information	46
1.1	Using this instruction for use	46
1.2	Icons	46
1.2.1	Device, packaging and accessories	46
1.2.2	Operating elements	47
1.3	Symbol convention	47
1.4	Glossary	48
1.5	Purpose	49
1.5.1	Intended use	49
1.5.2	Essential features	49
1.5.3	Applied parts	49
1.5.4	Indications	49
1.5.5	Contraindications	49
1.5.6	Restrictions on use	49
1.6	Basic safety instructions - ATTENTION!	50
1.7	User requirements	52
1.8	Information on product liability	52
1.9	Material compatibility	53
2	Product description	54
2.1	Whole view of the ME system	54
2.2	Product contents	54
2.3	Product properties	54
2.4	Information on the bacterial filter	55
2.5	Information on the battery	56
2.6	Warranty	56
3	Operation	57
3.1	Set-up and startup	58
3.1.1	Connecting the ASSKEA secretion aspirator proVisio® M28	58

3.2	Startup	59
3.2.1	Connecting the reusable secretion canister	59
3.3	Operation	62
3.3.1	Control elements and indicators of the ASSKEA secretion aspirator proVisio® M28	62
3.3.2	Setting the vacuum	64
3.3.3	Aspiration	64
3.4	Decommissioning	65
4	Maintenance	67
4.1	Cleaning and care in outpatient and inpatient care	67
4.1.1	General information	67
4.1.2	Cleaning and disinfection of the surface of the device	67
4.1.3	Cleaning and disinfection of the reusable secretion canister	68
4.1.4	Cleaning/disposal of the suction tube	68
4.1.5	Cleaning and disinfection of the rinsing bottle	68
4.2	Cleaning and care in the homecare area	69
4.2.1	General information	69
4.2.2	Cleaning of the surface of the device	69
4.2.3	Hygienic cleaning and disinfection of the reusable secretion canister	70
4.2.4	Hygienic cleaning and disinfection of the suction tube	70
4.2.5	Hygienic cleaning and disinfection of the rinsing bottle.....	71
4.3	Further use of the device	71
4.4	Maintenance and service	71
4.5	Inspection of the ASSKEA secretion aspirator proVisio® M28	72
5	Troubleshooting	73

6	Transport, storage and disposal	74
6.1	Decontamination prior to shipment	74
6.2	Storage.....	74
6.3	Disposal	74
7	Technical data	75
8	EMC information	76
8.1	Electromagnetic environment in which the ASSKEA secretion aspirator proVisio® M28 may be operated.....	76
8.2	Handling of electromagnetic interactions	79
8.3	List of all cables and transformers replaceable by the user	79
9	Order information accessories	80
10	Publishing information	81

1 User information

1.1 Using this instruction for use

Please read this entire instruction for use before operating the ASSKEA secretion aspirator proVisio® M28 for the first time. Alongside instructed medical professionals and trained relatives of patients, it is also intended that patients operate the device themselves, providing that they have also been instructed accordingly. In this case all functions can be used safely. If support is required for startup, operation or maintenance, please contact ASSKEA GmbH (see section 10). Please also report unexpected operating performance or incidents to ASSKEA GmbH.


Please read the safety instructions (see section 1.6) to avoid hazards.

This instruction for use forms an integral part of the ASSKEA secretion aspirator proVisio® M28. Therefore please keep this instruction for use in an easily accessible place.

Include this instruction for use when passing the ASSKEA secretion aspirator proVisio® M28 on to third parties..

1.2 Icons

1.2.1 Device, packaging and accessories

Symbol	Meaning
	Attention: possible bodily injury, health risks or possible property damage
	Follow the instruction for use
	Power supply unit
	Protection class II
	Protection class: Type BF (Body Floating)
	Temperature limitation
	Air pressure limitation
	Humidity limitation
	Do not use if packaging is damaged!
	International Protection: IP22 (see section 1.4)
	Note with useful information and tips

Symbol	Meaning
	This device must not be disposed of in domestic waste.
	Order number
	Serial number
	Lot number
	Date of manufacture
	Manufacturer
	Operating time (power-on time and power-off time)
	Do not reuse.

1.2.2 Operating elements

Symbol	Meaning
	Vacuum setting "Low" (vacuum \leq -0.2 bar)
	Vacuum setting "Medium" (-0.2 bar < vacuum \leq -0.5 bar)
	Vacuum setting "High" (-0.5 bar < vacuum \leq -0.9bar)
	LED status indicator for battery operation
	LED status indicator for mains operation
	On/off button

1.3 Symbol convention

Symbol	Meaning
•	Enumeration
1. 2.	Perform the actions in the specified order.

1.4 Glossary

A Aspirate Aspirate is the generic term for secretions, bodily fluids and flushing liquids that are typically accumulated during aspiration of the upper airways. These can be aspirated easily using the device described here.

C Contamination Contamination means that bacteria and viruses from the aspirate have come into contact with the device.

D Degree of protection The degree of protection indicates the protection of applied parts against electric shock. Applied parts of type BF have to be installed isolated from earth and are not suitable for direct application to the heart.

I incl. Abbreviation for "including"

IP22 International Protection/Protection Class
The Protection Class specifies the degree of protection of the device against contact and ingress of liquids.

The ASSKEA secretion aspirator proVisio® M28 is protected against finger access and falling water drops at an inclination of up to 15°.

M ME Medical electrical

MRI Abbreviation for "magnetic resonance imaging"
This technique generates sectional images of the human body with the aid of a very strong magnetic field for the purpose of analyzing the organs.

O Overflow Overflow means that aspirate is sucked into the interior of the device.

P Processing The processing procedure is required for each new patient. The term processing means that all parts that have (potentially) come into contact with the aspirate must be cleaned, disinfected and replaced if necessary.

Reconditioning may only be performed by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH.

R resp. Abbreviation for "respectively"

1.5 Purpose

1.5.1 Intended use

The ASSKEA secretion aspirator proVisio® M28 is a network-independent mobile medical aspirator and serves for temporary and primarily spontaneous aspiration of aspirate from the trachea.

Typical areas of application are:

- in outpatient and inpatient care (professional healthcare facility environment),
- in the homecare area (home healthcare environment), especially for the aspiration of aspirate from tracheostomy patients.

1.5.2 Essential features

- Generation of vacuum (high vacuum)
- Generation of flow (high flow)

1.5.3 Applied parts

There are no applied parts, however housing and suction tube are considered as applied parts of type BF.

1.5.4 Indications

- Tracheostomy patients
- Aspiration on impaired respiratory function
- Aspiration on blood, secretion and food residues from the oral cavity, the pharyngeal zone and the bronchial system

1.5.5 Contraindications

The ASSKEA secretion aspirator proVisio® M28 is contraindicated for the following applications:

- Liposuction
- Gynecology applications
- Dental applications
- Thoracic drainage
- Continuous drainage
- Applications in wound areas

1.5.6 Restrictions on use

- In medical rooms where potential equalization is required (e.g. heart surgery)
- In hazardous areas/in the MRI environment
- Outdoors/during transport

1.6 Basic safety instructions - ATTENTION!

Warning of damage due to improper power supply

Improper operation causes overvoltage in the device, which may be transmitted to the operator.

- Ensure prior to startup that the mains supply for connecting the ASSKEA secretion aspirator proVisio® M28 is suitable at 100 V to 240 V AC and a mains frequency of 50-60 Hz.
- Ensure prior to startup in UL listed markets such as the USA and Canada that the mains supply is designed for a supply voltage of 120 V AC.
- Only use the supplied power supply unit with the ASSKEA secretion aspirator proVisio® M28 (type: GTM91099-6015-3.0-T2 from GlobTek Inc.).

Risk of application under conditions that are not approved

- The device is not intended to be used in medical rooms where potential equalization is required (e.g. heart surgery).
- The device is not intended to be used in hazardous areas/in the MRI environment.
- The device is not intended to be used outdoors/during transport.

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never operate the device without bacterial filter.
- Prior to each application, make sure that the bacterial filter is dry and clean to ensure faultless operation. Pay attention that the bacterial filter does not get moist during the performance of cleaning/disinfection.
- Always perform the aspiration using a suitable sterile disposable catheter. The suction tube must never come into direct contact with the aspiration site.
- Follow the hygiene, cleaning and decontamination instructions.

Hazard of persons due to strangulation

- People may strangle themselves on the tubing or the mains cable, especially if tubes resp. cables are unduly long.
- Ensure that no unauthorized/uninvolved persons are near the device during aspiration.
- Store the device incl. accessories in the shipping carton until the next time of use.

Risk of damage due to electromagnetic phenomena

Medical electrical devices are subject to special precautionary measures regarding electromagnetic compatibility and must be installed and operated in accordance with the EMC information provided in the accompanying documents (see section 8).

Hazard of persons due to improper handling

- Only use the device for its intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. thoracic drainage).
- If aspiration procedures are carried out too often, slight bleeding may occur.
- When using the device for the treatment of children, the vacuum setting "Low" must be used.
- When using the power supply unit, ensure that it is connected to the aspirator firstly and afterwards to the mains supply (100 V to 240 V AC).
The disconnection of the power supply unit from the mains supply must be performed in exactly the opposite sequence (first remove the power supply unit from the mains supply (100 V to 240 V AC) and then remove it from the aspirator).
- Never touch parts of other devices than ME devices in the patient environment and the patient simultaneously.

Safety defects due to improper accessories and spare parts

The use of accessories and spare parts, detachable parts or materials that are not recommended by ASSKEA GmbH and specified in the instruction for use may compromise the safety and function of the device. Any warranty is excluded for damage caused by using accessories and spare parts not recommended or by improper use.

- Only use original accessories and spare parts.

Warning of safety defects due to improper connections of the ME system

The connection of the ME system with other devices or installations resp. pieces of equipment not recommended by ASSKEA GmbH and not specified in the instruction for use may compromise the safety and function of the ME system. Any warranty is excluded for damage caused by connecting devices or installations resp. pieces of equipment not recommended with the ASSKEA secretion aspirator proVisio® M28 or by improper use.

- Only connect recommended original parts with the ASSKEA secretion aspirator proVisio® M28.

Damage of the device by heat

- Do not cover the power supply unit.
- Keep the aspirator as well as the mains cable and the power supply unit away from other heat sources.
- Do not place the aspirator directly adjacent to other devices, since this can result in heavy warming of the aspirator.

Damage of the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Prior to each use check the housing for any damage and do not operate the device in case of apparent housing damage.
- Prior to each use check the suction tube, secretion canister and bacterial filter and if applicable any other accessories, that are subject to wear and damage, to ensure that the components are in a perfect condition and proper operation of the device is guaranteed. If this is not the case, replace the parts immediately.

Inspection of the internal power supply

Since the internal battery is not maintained automatically in a fully operational state, it is necessary to check the charging status regularly and if necessary, the battery must be replaced by service staff. The battery may only be replaced by trained service staff, since the replacement by inadequately trained personnel could result in a hazard (such as excessive temperatures, fire or explosion)!

Possible physiological effects and unobvious risks

- Select the vacuum range depending on the respective patient and on the medical indication to avoid personal injuries.
- Always place the device upright on a sturdy and flat base without slope. Make sure that the device cannot be knocked over or fall in such a way that persons could be hit by the falling device.
- Only operate the device with inserted vacuum pump and closed lid. The lid above the vacuum pump may only be opened by adequately trained service staff! Likewise, the vacuum pump may only be replaced by adequately trained service staff.
- Other devices, examinations or treatments may possibly be influenced by the device. For this reason, special attention should always be paid to other devices as well as to examinations and treatments performed in parallel so that any influence is detected as quickly as possible.
- When using the device, adequate lighting must always be ensured so that all labels can be identified clearly.
- Small, detached parts could be inhaled or ingested. Therefore ensure that no unauthorized persons, children or pets are near the device.
- Although the materials used have been tested for compatibility, in exceptional cases allergic reactions to the exposed materials on the device and its accessories may occur. This applies especially for contact injuries after prolonged contact. In such cases, turn to a doctor immediately.

Known, identifiable or foreseeable conditions for medical care within a domestic environment

- Children and pets should be kept away from the device. Make sure that the device cannot be knocked over or fall in such a way that persons could be hit by the falling device.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the domestic power supply.
- The device may not be used in damp rooms, bathrooms or showers. Avoid wet conditions on power supply unit, control panel and socket for the power supply unit.
- Ensure that fluff and dust are immediately removed from the device incl. accessories, in order to ensure unimpaired functionality. Furthermore no pests may be allowed in the proximity of the device, as there is a risk they may get inside the device causing damage.
- Never expose the device incl. accessories to direct sunlight as this could result in a strong heating effect, causing functional impairments.
- Some devices and sources that are commonly used in the domestic environment may cause potential faults in the device incl. accessories, e.g. fireplaces or radiant heaters (strong heating of the device) resp. inhalers or steam kettles (excessive air humidity). Do not operate such devices and sources in the vicinity of the ASSKEA secretion aspirator proVisio® M28.

Damage of the device by ingress of liquids

The ASSKEA secretion aspirator proVisio® M28 has the IP classification IP22 against the ingress of liquids. Nevertheless, protect the device from moisture..

- Do not use the device near splashing water.
- Do not use the device in damp rooms or while bathing or showering.
- Avoid wet conditions on power supply unit, control panel and socket for the power supply unit.
- Never submerge the device in water or other liquids (also not while it is switched off).

1.7 User requirements

The ASSKEA secretion aspirator proVisio® M28 may only be operated and used by instructed and appropriately trained persons. Familiarize yourself with the mode of operation of the ASSKEA secretion aspirator proVisio® M28 prior to startup.

Training on how to handle the ASSKEA secretion aspirator proVisio® M28 is provided by ASSKEA GmbH or an authorized distribution partner of ASSKEA GmbH. Product training takes approximately one to two hours and, among other things, includes an explanation of the design and function of the device, the handling of the device, the cleaning and disinfection as well as the procedure to be followed for each new patient and for disposal.

The training should be repeated regularly every 24 months. Each participant receives a certificate as proof of training.

1.8 Information on product liability

The liability for the operation of the device is channeled to the operator in the following cases:

- the device is used outside its intended use,
- the device is not used in accordance with the instruction for use,
- the device is opened by unauthorized personnel,
- the safety seal is removed/broken,
- installation, settings, enhancements, routine maintenance or repairs are performed by unauthorized personnel,
- original accessories and spare parts have not been used.

Advice for the responsible organization:

The assembly of ME systems and modifications during their actual service life shall be evaluated based on the requirements of the applying standards.

1.9 Material compatibility

Aggressive substances can damage the device and the accessories.

Follow the instructions for cleaning and care (see section 4.1 resp. 4.2).

2 Product description

2.1 Whole view of the ME system


Fig. 1

- A Secretion canister (reusable) with integrated overflow protection in the lid (fig. 19)
- B Control panel
- C On/off button
- D Socket for power supply unit (12V DC)
- E Power supply unit (GlobTek GTM91099-6015-3.0-T2)
- F Mains cable
- G Suction tube
- H Lid vacuum pump
- I Carrying handle

2.2 Product contents

- the ASSKEA secretion aspirator proVisio® M28
- the instruction for use of the ASSKEA secretion aspirator proVisio® M28
- the power supply unit GlobTek GTM91099-6015-3.0-T2
- reusable secretion canister system (incl. suction tube and fingertip)
- rinsing bottle, reusable (250 ml)
- charging instructions, multilingual
- safety instruction for handling battery packs
- "Used Medical Device" label and decontamination certificate
- test report in accordance with IEC 62353
- optional accessories (depending on the order)

2.3 Product properties

Hazards of persons due to improper handling

- Only use the device for its intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. thoracic drainage).

Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Do not operate the device in case of apparent housing damage.

The ASSKEA secretion aspirator proVisio® M28 is an especially handy, battery powered medical aspiration device for mobile use.

The ASSKEA secretion aspirator proVisio® M28 has a flow rate of max. 30 l/min ± 2 l/min (see section 7 "Technical data").

The ASSKEA secretion aspirator proVisio® M28 is operated via the internal battery or via the supplied power supply unit that can also be used to recharge the battery. A rotary vane pump generates the desired vacuum.

After 5 minutes of operation, the device switches off automatically. It is recommended to leave the device switched off for 5 minutes; however, the device can be switched on again immediately, as needed.

Moreover an overheating protection prevents the battery from overheating by interrupting the charging process if the temperature of the battery is too high (e.g. due to unfavorable environmental conditions).

After it is switched on, the vacuum pump generates a vacuum in the secretion canister, which is used to aspirate secretion (via a suitable sterile aspiration set). The aspirate is directed away from the patient and collected in the secretion canister. If the secretion canister is full, the aspiration is stopped by means of a mechanical overflow protection integrated into the lid of the secretion canister.

The desired vacuum is set by buttons in the control panel (fig. 1(B), see section 1.2.2)). The secretion canister is attached by sliding it directly into the housing structure of the aspirator (fig. 8, 9). This also establishes the vacuum connection between the aspirator and secretion canister.

The device is charged or operated via the supplied 12V power supply unit (fig. 1(E)).

2.4 Information on the bacterial filter

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic microbes in the suction material can pose a risk to health.

- Never operate the device without bacterial filter.
- Prior to each application, make sure that the bacterial filter is dry and clean to ensure faultless operation.

The bacterial filter integrated between the vacuum pump and the secretion canister consists of a hydrophilic (self-sealing) material and has an efficiency of 99% against bacteria and viruses. Thus it effectively protects the vacuum pump against overflow and contamination and enables fast, simple and cost-effective processing.

The **integrated bacterial filter** is only intended for **use on one patient**. It must therefore be replaced by the operator each time a new patient uses the device (see section 3.2.1). If the device is used on one patient exclusively, the bacterial filter must be replaced at least every 2 weeks.

To ensure consistent performance, the integrated bacterial filter must also be replaced by authorized service staff if it has come into **contact with the aspirate** (blocked) or **during maintenance/repair**.

Since the drive unit is mechanically isolated from the vacuum pump, it cannot become contaminated during use independent from the bacterial filter.

2.5 Information on the battery

The charge state of the battery is indicated by the corresponding LED status indicator (for explanation see section 3.3.1).

Prior to the initial startup of the ASSKEA secretion aspirator proVisio® M28, it is urgently recommended to charge the battery completely and to repeat this after the first uses.

The ASSKEA secretion aspirator proVisio® M28 is equipped with a lithium ion battery which only has a limited self-discharge on its disposal in contrast to conventional battery types.

Lithium ion batteries do not have a memory effect. Therefore they can and should be recharged at any time after initial charging. Only frequent short charging should be avoided. The typical service life of the battery amounts to approximately 500 load cycles, subsequently approximately 70% of the primary capacity are still available.

The battery of the ASSKEA secretion aspirator proVisio® M28 is protected against depth discharge by means of protective measures; however the instructions on charging mentioned above should be observed. The battery is also protected against overheating during charging. If the battery temperature is exceeded during charging due to ambient conditions not intended, charging will be automatically interrupted temporarily to allow cooling. This feature ensures the safety of the battery and protects it from wear.

2.6 Warranty

WARNING:
Modifications of the ME device are not permitted.

The ASSKEA secretion aspirator proVisio® M28 is covered by warranty for 2 years. No extension or renewal of the warranty is granted in the event that warranty work is performed.

The battery is covered by warranty for 6 months. Wearing parts are excluded from the warranty.

ASSKEA GmbH is responsible for impacts on safety, reliability and specific performance only if:

- original ASSKEA accessories and replacement parts are used,
- maintenance and repair are performed by professionals authorized by ASSKEA GmbH or by ASSKEA GmbH itself,
- the ASSKEA secretion aspirator proVisio® M28 is used and operated in accordance with the instruction for use and within its intended use.

All warranty claims shall be void if:

- the device is opened by unauthorized personnel,
 - the safety seal is removed/broken,
 - repairs are performed by unauthorized personnel,
 - modifications are made to the device,
- since the basic safety of the device can no longer be guaranteed in these cases and functional limitations may occur.

3 Operation

Hazard of persons due to improper handling

- Please read through sections 3.1 and 3.2 carefully!
- Only perform aspiration in the respiratory tract after instruction by specialized personnel!
- Only use suitable sterile disposable catheters for aspiration! Commercially available sterile disposable catheters from CH10 to CH16 that are intended for the use in the tracheal area can be connected.
- If aspiration procedures are carried out too often, slight bleeding may occur.
- When using the device for the treatment of children, the vacuum setting "Low" must be used.

Malfunction due to aspirated secretion

- Change the reusable secretion canister when it is half-full to curtail frothing. If the reusable secretion canister is full, the mechanical overflow protection is triggered. This disrupts the aspiration process.
- Switch off the device when emptying the reusable secretion canister.
- If secretion has been aspirated into the device, it must be properly processed by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH!

Damage of the device due to insufficient acclimatization

- After the device has been exposed to temperatures between -25°C and $+70^{\circ}\text{C}$ during transport/storage according to the technical data (see section 7), it must firstly acclimatize for ca. 2 h at room temperature (ca. 20°C) before the intended use is possible.

3.1 Set-up and startup

The following section describes the operating elements, connections and the startup of the ASSKEA secretion aspirator proVisio® M28:

3.1.1 Connecting the ASSKEA secretion aspirator proVisio® M28

Known or identifiable conditions for medical care within a domestic environment

- Children and pets should be kept away from the device to ensure that it cannot be knocked over or dropped.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the domestic power supply.
- The device may not be used in damp rooms, bathrooms or showers. Avoid wet conditions on power supply unit, control panel and socket for the power supply unit. Never submerge the device in water or other liquids (also not while it is switched off).

Check the power supply unit incl. mains cable for potential damages before each use and replace it if there is any damage.

Use the socket for the 12V power supply unit (fig. 1(D)) to connect the device to the mains power supply with the aid of the supplied power supply unit to charge or operate it as needed. Make sure that the device is positioned in such a way that it is possible to easily disconnect it later.

Information on the permissible environmental conditions during operation can be found in section 7 "Technical data".


Fig. 2



Connection of the ASSKEA secretion aspirator proVisio® M28 to patient and accessories

3.2 Startup

- Take the device and the accessories out of the packaging.
- Always place the device on a sturdy and flat surface without slope.
- Prior to the initial startup it is essential to follow the safety instructions in section 1.6.
- Only operate the device with inserted bacterial filter and connected secretion canister. Bacterial filter and secretion canister are already fitted on delivery. Note the following instructions for all further applications.
- Always have one backup bacterial filter at hand, since it is absolutely essential for safe operation!
- Prior to each application, check that the reusable secretion canister is completely clean to avoid frothing.


3.2.1 Connecting the reusable secretion canister

1. Fold the moving part of the secretion canister lid (visioClick®) slightly upwards and place the secretion canister lid on the secretion canister.



3 Operation


2. Fold the moving part of the secretion canister lid (visioClick®) downwards so that the hooks snap into place and the secretion canister is securely closed.

3. Insert the secretion canister into the canister fitting provided on the device by guiding it over the rack on the bottom and pivoting it into the device until it snaps into place.



4. Connect the suction tube to the port in the middle of the secretion canister lid. In doing so, ensure that the tube is fixed sufficiently tight.

5. Connect the other ending of the suction tube to the aspiration catheter with the aid of a fingertip.

Insertion of the bacterial filter


Insert the bacterial filter into the device by introducing the filter with the lever pointing vertically downwards in the opening provided. Subsequently move the lever a quarter turn clockwise so that the lever is in the horizontal position afterwards.

Please note that the reusable secretion canister cannot be connected if the bacterial filter has not been inserted correctly.

3.3 Operation

3.3.1 Control elements and indicators of the ASSKEA secretion aspirator proVisio® M28

- A Vacuum setting "High"
- B Vacuum setting "Medium"
- C Vacuum setting "Low"
- D LED status indicator for battery operation
- E LED status indicator for mains operation
- F On/off button

Fig. 11

Symbol	Description	Meaning of display
	button illuminates continuously green	device is running, vacuum setting "High" is adjusted (-0.5 bar < vacuum ≤ -0.9 bar)
	button illuminates continuously red	canister full/filter blocked (motor switches off automatically)
	button illuminates continuously green	device is running, vacuum setting "Medium" is adjusted (-0.2 bar < vacuum ≤ -0.5 bar)
	button illuminates continuously red	canister full/filter blocked (motor switches off automatically)
	button illuminates continuously green	device is running, vacuum setting "Low" is adjusted (vacuum ≤ -0.2 bar)
	button illuminates continuously red	canister full/filter blocked (motor switches off automatically)
	display illuminates continuously green	battery capacity 100 % to 60 %
	display illuminates continuously orange	battery capacity < 60 % to 10 %
	display illuminates continuously red	Battery capacity < 10%, connect power supply unit immediately. Switch off the device after completion of the aspiration procedure for charging!
	display blinks green	Battery is being charged.
	display blinks red	Overtemperature protection, interruption of charging process until battery temperature returns to permissible range.
		Battery defective, operation only possible with connected power supply unit. Please contact your service partner!
	display illuminates continuously green	Power supply unit is connected, voltage available.
	display blinks red	Incorrect power supply unit is connected, no charging of battery possible!
	button illuminates continuously green	Device is switched on and ready for use.
	button illuminates continuously orange	Device still runs for 30 sec. before it switches off automatically. Afterwards the button will still be illuminated for 10 seconds.
	button illuminates continuously red	Internal error, operation is not possible. Please contact your service partner!

3 Operation

3.3.2 Setting the vacuum

1. Switch on the device by pressing the on/off button (fig. 11(F)).
2. Select the desired vacuum range by pressing the corresponding button ("Low", "Medium" or "High") on the control panel (fig. 11 (A, B, C)).

As described in section 3.3.1, the vacuum setting buttons correspond to certain vacuum ranges.

3. You can perform fine adjustments of the set vacuum via the additional air opening on the fingertip (fig. 12).


Fig. 12

4. Switch off the device by holding the on/off button (fig. 11(F)) pressed for 1-2 seconds.

3.3.3 Aspiration

1. Perform the aspiration procedure as you were instructed by specialized personnel.
2. Control the aspiration procedure and the desired vacuum via the additional air opening on the fingertip (fig. 12).
3. Decommission the device after the aspiration procedure according to section 3.4.

In any case, please take note of section 2.4 and the important information on the handling of the bacterial filter! Ensure that the battery has a sufficient capacity during the aspiration procedure. The battery status display must signalize green or orange. If the display turns red, the power supply unit must be plugged in immediately. Only in this case the aspiration can be continued without any interruption.

3.4 Decommissioning

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Replace the bacterial filter at least every 2 weeks when used exclusively for one patient!
- Wear suitable disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- For each new patient, processing by ASSKEA GmbH or by an authorized service partner of ASSKEA GmbH is strictly required!
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The aspirate and the parts contaminated with aspirate must be disposed of properly.

1. Switch off the device after the aspiration by holding the on/off button (fig. 11(F)) pressed for 1-2 seconds.
2. Disconnect the power supply unit from the mains supply (100 V to 240 V AC) if attached and remove the device plug from the aspiration device afterwards.
3. Remove the suction tube (fig. 1(G)) from the secretion canister.
4. Clean or dispose of the suction tube properly (see section 4.1.4 resp. 4.2.4).

5. Remove the reusable secretion canister from the device by pressing the unlocking mechanism on the secretion canister lid and pivoting the secretion canister out of the device afterwards.


Fig. 13


Fig. 14

3 Operation

6. Fold the moving part of the secretion canister lid (visioClick®) upwards so that the hooks are released and take off the lid from the secretion canister.

Fig. 15

7. Empty and clean the reusable secretion canister and the lid according to section 4.1.3 resp. 4.2.3.
8. Remove the used bacterial filter before a new patient uses the device or when it has already been used for one patient for 2 weeks. Move the lever that is in the horizontal position a quarter turn counterclockwise so that the lever points vertically downwards subsequently.


Fig. 16


Fig. 17

9. When the lever is in this position, pull the bacterial filter out of the fitting and dispose of the bacterial filter properly.


Fig. 18

10. Insert a new bacterial filter (see section 3.2.1).
11. Clean the surface of the device according to section 4.1.2 resp. 4.2.2.
12. Store the device in the shipping carton until next operation.

4 Maintenance

4.1 Cleaning and care in outpatient and inpatient care

4.1.1 General information

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never operate the device without bacterial filter.
- Replace the bacterial filter at least every 2 weeks when used exclusively for one patient!
- Wear suitable disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- For each new patient, processing by ASSKEA GmbH or by an authorized service partner of ASSKEA GmbH is strictly required!
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The aspirate and the parts contaminated with aspirate must be disposed of properly.

Health risks due to the handling of disinfectants

- The use of appropriate protective clothing during disinfection is recommended.
- Follow the disinfectant manufacturer's instructions.

Damage of the device due to improper cleaning agents

- Do not use disinfectants containing acetone. These may damage or disfigure the housing parts as well as the secretion canister lid.
- Follow the instruction for use provided by the manufacturers of the used disinfectants, particularly concerning the compatibility with materials and surfaces as well as the concentration information.
- ASSKEA GmbH recommends "Sekusept® aktiv" for immersion disinfection of the accessories and "Incidin® Plus" or "Incidin® Liquid" for wipe disinfection of the device.

4.1.2 Cleaning and disinfection of the surface of the device

Clean the surfaces of the device regularly and disinfect them daily.

- The device can be wiped with a damp, lint-free cloth.
- Follow the instructions in section 4.1.1 for wipe disinfection.

Repeated cleaning and disinfection procedures may result in minor discolorations of the plastic components of the housing. However these do not affect the function of the device.

If the interior of the device comes into direct contact with liquids or solids, the device must be inspected by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH.

4.1.3 Cleaning and disinfection of the reusable secretion canister


Fig. 19

- A Reusable secretion canister lid
B Overflow protection with sphere

1. Empty the reusable secretion canister and dispose of the aspirate properly.
2. Remove the sphere before performing the following cleaning instructions (fig. 19(B)).
3. Rinse all components of the reusable secretion canister under running water.
4. Immerse all components of the reusable secretion canister in disinfectant solution, considering the specific concentration quantity.
5. Afterwards, rinse the components thoroughly and let them dry.
6. The reusable secretion canister can also be autoclaved for 20 minutes at 121°C, in that case please remove the overflow protection in advance (fig. 19(B)).
7. After autoclaving, fix the overflow protection in the lid considering the orientation of the opening (fig. 19(B)).

8. Reinsert the sphere according to fig. 19.

ASSKEA GmbH recommends to replace the canister system **every 4 weeks** if it is used and disinfected frequently!

4.1.4 Cleaning/disposal of the suction tube

1. Dispose of the suction tube properly before the device is used by a new patient!
2. If the device is used for one patient only, rinse the suction tube with clear water using the rinsing bottle after each aspiration and place it daily in the disinfectant solution recommended by ASSKEA GmbH for immersion disinfection (see section 4.1.1).
3. Afterwards, rinse the suction tube thoroughly with clear water and let it dry.

Even when the suction tube is used for one patient only, ASSKEA GmbH recommends to replace the suction tube at least **every 4 weeks** since the material may become discolored and brittle due to frequent applications!

4.1.5 Cleaning and disinfection of the rinsing bottle

1. Remove the lid of the rinsing bottle and rinse both parts under running water.
2. Immerse the rinsing bottle and the lid in disinfectant solution, considering the specific concentration quantity.
3. Afterwards, rinse the components thoroughly and let them dry.
4. Place the lid back on the rinsing bottle.

ASSKEA GmbH recommends to replace the rinsing bottle **every 4 weeks** if it is used and disinfected frequently!

4.2 Cleaning and care in the homecare area

4.2.1 General information

Health risks due to the handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate cause health risks.

- Never operate the device without bacterial filter.
- Replace the bacterial filter at least every 2 weeks when used exclusively for one patient!
- Wear suitable disposable gloves in doing so.
- The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- Prior to passing the device on to another patient resp. to other persons, processing by ASSKEA GmbH or by an authorized service partner of ASSKEA GmbH is strictly required! For that purpose observe the information regarding the further use of the device (see section 4.3).
- After each aspiration, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- The aspirate and the parts contaminated with aspirate must be disposed of properly.

Health risks resp. damages of the device due to the handling of disinfectants

- If you have the possibility to use disinfectants, please observe the information in section 4.1. Otherwise please follow the indications of section 4.2.2 till 4.2.5!

4.2.2 Cleaning of the surface of the device

Clean the surfaces of the device regularly.

The device can be wiped with a damp, lint-free cloth. Repeated cleaning procedures may result in minor discolorations of the plastic components of the housing. However these do not affect the function of the device.

If the interior of the device comes into direct contact with liquids or solids, the device must be inspected by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH.

4.2.3 Hygienic cleaning and disinfection of the reusable secretion canister


Fig. 20

- A Reusable secretion canister lid
- B Overflow protection with sphere

1. Empty the reusable secretion canister and dispose of the aspirate properly.
2. Remove the sphere before performing the following cleaning instructions (fig. 20(B)).
3. Rinse all components of the reusable secretion canister under running water.
4. Place all components of the reusable secretion canister in a water bath at 65°C temperature for 15 minutes once a day.
5. Afterwards, rinse the components thoroughly and let them dry.
6. Reinsert the sphere according to fig. 20.

ASSKEA GmbH recommends to replace the canister system **every 4 weeks** if it is used and cleaned/disinfected hygienically frequently!

4.2.4 Hygienic cleaning and disinfection of the suction tube

1. Dispose of the suction tube properly before the device is used by a new patient!
2. If the device is used for one patient only, rinse the suction tube with clear water using the rinsing bottle after each aspiration and place it in a water bath at 65°C temperature for 15 minutes once a day.
3. Afterwards, rinse the suction tube thoroughly with clear water and let it dry.

Even when the suction tube is used for one patient only, ASSKEA GmbH recommends to replace the suction tube at least **every 4 weeks** since the material may become discolored and brittle due to frequent applications!

4.2.5 Hygienic cleaning and disinfection of the rinsing bottle

1. Remove the lid of the rinsing bottle and rinse both parts under running water.
2. Place the rinsing bottle and the lid in a water bath at 65°C temperature for 15 minutes once a day.
3. Afterwards, rinse the components thoroughly and let them dry.
4. Place the lid back on the rinsing bottle.

ASSKEA GmbH recommends to replace the rinsing bottle **every 4 weeks** if it is used and cleaned/disinfected hygienically frequently!

4.3 Further use of the device

The ASSKEA secretion aspirator proVisio® M28 is suitable for further use. However, prior to passing the device on to other patients resp. persons, it has to be processed professionally. For this purpose, please hand the ASSKEA secretion aspirator proVisio® M28 over to ASSKEA GmbH or to qualified personnel authorized by ASSKEA GmbH. In this regard, please observe the indications in section 6.1!

4.4 Maintenance and service

Risk of impairments of performance due to aging effects

Only use the device within its indicated service life of 5 years. If you notice impaired performance or other problems, follow the instructions on troubleshooting in section 5 and contact the service if necessary.

Risk due to maintenance and service during operation

Maintenance and service must not be performed while the device is in use! The device must be switched off before maintenance and service work are performed.

Risk due to the performance of maintenance and service by unauthorized persons

Maintenance and service may only be performed by specialized personnel authorized by ASSKEA GmbH or by ASSKEA GmbH itself.

Risk due to the performance of maintenance and service without sufficient documents

Maintenance and service may only be performed considering the service instructions for the ASSKEA secretion aspirator proVisio® M28. ASSKEA GmbH will also provide additional documents on request, in case these can support the service staff.

WARNING: Modifications of the ME device are not permitted.

- The ASSKEA secretion aspirator proVisio® M28 is maintenance-free if used in accordance with the instruction for use, exempt from components with limited durability.
- Perform a visual and functional check prior to each application. Also include the accessories of the device in doing so.
- The device may only be opened and repaired by ASSKEA GmbH or by specialized personnel authorized by ASSKEA GmbH in compliance with the service documents specified by the manufacturer as well as with technical and hygienic safety measures.
- The device may be returned to ASSKEA GmbH for repair either directly or via the specialty dealer you purchased the device from.
- All accessories must be cleaned and disinfected before the device is returned. The device itself may be treated with a surface disinfectant.
- In order to avoid time lags, notify ASSKEA GmbH prior to returning the device and mark the shipping carton of the device with the label "Used Medical Device".

The corresponding product return form with the "Used Medical Device" label can be found on the web page of ASSKEA GmbH at www.asskea.de in the service area under product return.

ASSKEA GmbH neither ensures proper operation for the ASSKEA secretion aspirator proVisio® M28 nor is ASSKEA GmbH liable for any resulting property damage or personal injuries if:

- original ASSKEA accessories and spare parts are not used,
- the user instructions in this instruction for use are not followed,
- installation, settings, modifications, enhancements or repairs are not performed by ASSKEA GmbH or by personnel authorized by ASSKEA GmbH.

4.5 Inspection of the ASSKEA secretion aspirator proVisio® M28

ASSKEA GmbH offers its partners and customers fast and professional processing as well as required testing services.

5 Troubleshooting

Malfunction	Display	Probable cause	Remedy
Device does not start	  do not illuminate	Battery is discharged.	Connect the power supply unit for charging.
	 blinks red	Battery is discharged and impermissible input voltage.	Connect the correct power supply unit, check the power supply.
	 blinks red	Internal error	Please contact service!
Device does not charge	 blinks red and  illuminates green	<ul style="list-style-type: none"> Battery temperature too high for charging Battery is defective. 	<ul style="list-style-type: none"> No measure required, automatic resumption of charging when battery has cooled down. Please contact service!
	 blinks red	Impermissible input voltage	Connect the correct power supply unit, check the power supply.
	 illuminates green/orange and/or  illuminates green	<ul style="list-style-type: none"> Bacterial filter is blocked. Leak in the suction line or in the canister Service life of vacuum pump expired 	<ul style="list-style-type: none"> Replace the bacterial filter. Check that secretion canister lid and tube connections are tight and no leakages are apparent. Please contact service!
No flow rate	   continuous red light	<ul style="list-style-type: none"> Float sphere of the overflow protection blocks the aspiration lid inlet. Bacterial filter is blocked. Suction tube is blocked. Bacterial filter is not inserted. 	<ul style="list-style-type: none"> Check fluid level in the secretion canister, empty the canister. Replace the bacterial filter. Replace the suction tube. Insert the bacterial filter.
Device switches off	 blinks red	Internal error	Please contact service!

6 Transport, storage and disposal

6.1 Decontamination prior to shipment

Before the ASSKEA secretion aspirator proVisio® M28 is passed on, the device must be properly processed by ASSKEA GmbH or by specialized personnel authorized by ASSKEA GmbH to protect subsequent users. Processing is mandatory in accordance with the German Medical Device Operator Ordinance (MPBetreibV), German Medical Devices Act (MPG) and the manufacturer's instructions.

ASSKEA GmbH offers its partners and customers fast and professional processing as well as required testing services.

Before the ASSKEA secretion aspirator proVisio® M28 is shipped to ASSKEA GmbH, it must be cleaned and disinfected. Please follow the instructions in section 4.1 resp. 4.2! Please affix the supplied "Used Medical Device" label to the shipping carton! Please notify ASSKEA GmbH of the return in advance. The product return form can be found on the web page of ASSKEA GmbH at www.asskea.de in the service area under product return.

6.2 Storage

Store the ASSKEA secretion aspirator proVisio® M28 in accordance with the indications in the technical data (section 7)!

Charge the battery of the ASSKEA secretion aspirator proVisio® M28 before storing the device. This ensures that it is operational at all times.

If the ASSKEA secretion aspirator proVisio® M28 is not used for a longer period of time (approximately 10 months), the battery must be fully recharged again!

6.3 Disposal

- The components of the device must be disposed of properly at the end of the product's service life.
- Ensure that the disposed components are clean and carefully sorted by material.
- The housing material is labeled with a material symbol mark and is fully recyclable.
- Decontaminate the device and the accessories prior to disposal.
- According to EU Directives 2012/19/EU Waste Electrical and Electronic Equipment (WEEE) and 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II), the device must not be disposed of in domestic waste.
- The device and the accessories may be disposed of via ASSKEA GmbH or the service partner.
- In relation to proper disposal of potential biohazards, contact your local authority.
- Outside the EU: Follow the disposal requirements of your country!

7 Technical data

Model name	ASSKEA proVisio® M28
Flow rate* (measuring point at suction tube nozzle)	30 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg) (high vacuum) Conversion factor: 1 kPa ~ 7.5 mmHg; 1 kPa ~ 10 mbar
Canister	Reusable secretion canister (1000 ml)
Suction tube (item number: 100713-2)	PVC, inner diameter 6 mm, wall thickness 2 mm, length 150 cm, vacuum stable up to -0.9 bar
Nominal voltage of the power supply unit	In: AC 100 – 240 V ~ 50-60 Hz/1.5 A Out: DC 12 V/5.0 A
Maximum load current	5.0 A
Permissible input voltage	12 V
Power consumption at 12 V	60 W
Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pursuant to 93/42/EEC, IX	Ila
Protection class pursuant to IEC 60601-1	II
IP degree of protection	IP22
CE marking	CE1434
Sound emission	without vacuum: approximately 60 dB(A)

Ambient conditions

Transport/Storage

ambient temperature: -25 °C to +70 °C
relative humidity:
up to 90 %, non-condensing
air pressure: 700 hPa to 1060 hPa

Operation

ambient temperature: +5 °C to +40 °C
relative humidity:
15 % to 90 %, non-condensing
air pressure: 700 hPa to 1060 hPa

Battery	10.8 V; 5800 mAh (lithium ion battery)
Charging time if battery is empty	~ 3 h
Charging time if battery is approximately 50% full	~ 2 h
Power supply unit	GTM91099-6015-3.0-T2, GlobTek Inc.
Dimensions basic device (H x W x D)	213 mm x 260 mm x 124 mm
Weight (basic device)	2.8 kg
Operating time	5 minutes ON; 5 minutes OFF
Runtime in battery operation	approximately 45 minutes, depending on the runtime of the motor
Lifetime	5 years
Item number (REF)	100739

* The information provided may differ depending on the altitude above sea level, the prevailing air pressure and the air temperature.

8 EMC information

WARNING: The use of the ASSKEA secretion aspirator proVisio® M28 directly adjacent to or stacked with other devices should be avoided, since this could lead to impermissible operation. If the use of the ASSKEA secretion aspirator proVisio® M28 adjacent to or stacked with other devices is required, the ASSKEA secretion aspirator proVisio® M28 and the other devices should be monitored in order to verify proper operation in this arrangement!

WARNING: The use of accessories and spare parts, transformers and cables for the ASSKEA secretion aspirator proVisio® M28 not indicated or provided by ASSKEA GmbH may increase the electromagnetic emissions or reduce the electromagnetic immunity of the ASSKEA secretion aspirator proVisio® M28, resulting in impermissible operation. No warranty is provided for damages caused by using accessories and spare parts, transformers and cables not recommended or by improper use. Only use original ASSKEA accessories and spare parts!

WARNING: The use of the indicated or provided accessories and spare parts, transformers and cables for other devices than the ASSKEA secretion aspirator proVisio® M28 may increase the electromagnetic emissions or reduce the electromagnetic immunity. No warranty is provided for damages caused by using the indicated or provided accessories and spare parts, transformers and cables with other devices or by improper use. Only use the accessories and spare parts, transformers and cables with the ASSKEA secretion aspirator proVisio® M28!

WARNING: Portable and mobile RF communication equipment (incl. peripheral devices such as antenna cables and external antennas) may influence medical electrical devices and therefore should not be used within a range of 30 cm of any part of the ASSKEA secretion aspirator proVisio® M28 incl. its cables. Otherwise the performance of the device may be impaired.

WARNING: The ASSKEA secretion aspirator proVisio® M28 may influence other devices, examinations and treatments electromagnetically. For this reason, special attention should always be paid to other devices and examinations or treatments performed in parallel so that any influence is detected as quickly as possible.

The ASSKEA secretion aspirator proVisio® M28 meets the requirements of IEC 60601-1-2/EN 60601-1-2 "Electromagnetic Compatibility - Medical Electrical Equipment" without deviations or restrictions. Therefore electromagnetic interference is reduced to a minimum. Follow the indicated instructions and guidelines to sustain the basic safety and the essential features of the ASSKEA secretion aspirator proVisio® M28 over its entire service life.

8.1 Electromagnetic environment in which the ASSKEA secretion aspirator proVisio® M28 may be operated

The ASSKEA secretion aspirator proVisio® M28 is intended for operation in the electromagnetic environment specified below, in which the RF disturbances are controlled. The customer or the user of the ASSKEA secretion aspirator proVisio® M28 must ensure that it is operated in such an environment.

The environments for intended operation include professional health-care institutions and the home care environment. Special environments, such as the proximity of RF surgery or MRI resp. environments in which the intensity of the EMC disturbances is high, are excluded.

Emission limits

Conducted and radiated RF emissions	CISPR 11
-------------------------------------	----------

Harmonic distortion	IEC 61000-3-2
---------------------	---------------

Voltage fluctuations/flicker	IEC 61000-3-3
------------------------------	---------------

Housing

Phenomenon	Test method	Immunity test levels
------------	-------------	----------------------

Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air
-------------------------	---------------	-----------------------------

Radiated RF disturbances	IEC 61000-4-3	28 V/m 80 MHz to 6 GHz 80% AM at 1 kHz
--------------------------	---------------	--

Magnetic field at rated power frequency (50 Hz)	IEC 61000-4-8	30 A/m
---	---------------	--------

Power supply AC

Phenomenon	Test method	Immunity test levels
------------	-------------	----------------------

Fast transient electrical disturbances/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
---	---------------	---------------------------------------

Impulse voltages/surges	IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground
-------------------------	---------------	--

Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz
---------------------------	---------------	---------------------------

		6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
--	--	---

Voltage dips	IEC 61000-4-11	0% U_T for ½ period At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T for 1 period 70% U_T for 25 periods one-phase at 0°
--------------	----------------	---

Voltage interruptions	IEC 61000-4-11	0% U_T for 250 cycles
-----------------------	----------------	-------------------------

Note: U_T is the AC mains voltage prior to the application of the test levels.

Power supply DC

Phenomenon	Test method	Immunity test levels
Fast transient electrical disturbances/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Impulse voltages/surges	IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Electrical transient conduction along supply lines	ISO 7637-2	see ISO 7637-2

Patient connection

Phenomenon	Test method	Immunity test levels
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

Signal input and output channels

Phenomenon	Test method	Immunity test levels
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air
Fast transient electrical disturbances/bursts	IEC 61000-4-4	±1 kV 100 kHz repetition frequency
Impulse voltages/surges	IEC 61000-4-5	±2 kV line-to-ground
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

8.2 Handling of electromagnetic interactions

Although electromagnetic interference of the ASSKEA secretion aspirator proVisio® M28 has been reduced to a minimum, electromagnetic disturbances between the ASSKEA secretion aspirator proVisio® M28 and other devices cannot be excluded. Therefore you should comply with the specified requirements and instructions regarding the permissible electromagnetic environment at any case and monitor the ASSKEA secretion aspirator proVisio® M28 in order to ensure proper operation and to prevent adverse events for patients and users. Select another location for the ASSKEA secretion aspirator proVisio® M28, if the permissible electromagnetic environment cannot be ensured or if you have observed functioning of the ASSKEA secretion aspirator proVisio® M28 or of another device in the vicinity that is not intended.

Since electromagnetic propagation is affected by absorption and reflection from structures, objects and people and since the field strength from fixed transmitters cannot be predicted in advance with accuracy, a site survey of the electromagnetic phenomena at the location at which the ASSKEA secretion aspirator proVisio® M28 shall be operated should be considered to assess the existing electromagnetic environment at this location. If unusual performance is detected, additional measures may be required, for example a modified orientation or another location of the ASSKEA secretion aspirator proVisio® M28. If the essential features of the ASSKEA secretion aspirator proVisio® M28 are impaired by electromagnetic disturbances, a reduced or non-existent flow rate and vacuum may be expected.

8.3 List of all cables and transformers replaceable by the user

Name	Specification	Maximum length
Power supply unit incl. connecting cable	Type: GTM91099-6015-3.0-T2 Manufacturer: GlobTek Technical data: 100-240 V~, 50-60 Hz, 1.5 A (in) 12 V DC, 5 A (out)	1.20 m
Mains cable	Type: H03VVH2-F Manufacturer: HAWA Technical data: 250 V~, 2.5 A	1.80 m

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