ASSKEA proVisio[®]Basic
ASSKEA proVisio[®]Basic battery

REF 100738-1

REF 100738-2

de

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en

Instruction for Use

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Table of Contents

1	Important Information	48
2	Intended Purpose	49
3	Safety Instructions	53
4	Warranty and Product Liability	56
5	Symbols and Terms	57
5.1	Symbols on Device, Packaging and Accessories	57
5.2	Glossary	58
6	Product Description	59
6.1	Device and Accessories - Overview	59
6.2	Scope of Delivery	61
6.3	Applied Parts	61
6.4	Connection between Device and Accessories	62
6.5	Battery	63
6.6	Reusable Secretion Canister System	63
6.7	Disposable Secretion Canister System	64
7	Start-up	65
7.1	Starting up the device	65
7.2	Connecting the Reusable Secretion Canister System	65
7.3	Connection of the Disposable Secretion Canister System	67
8	Operation	69
8.1	Control and Display Elements	69
8.2	Operation and Function	70
8.3	Functional Check	70
8.4	Performing Suction	71
9	Follow-up	71
9.1	Reusable secretion canister system is full and must be emptied	71
9.2	Disposable secretion canister system is full and must be replaced	72
10	Cleaning and Disinfection	73
10.1	WARNINGS	73
10.2	Device Surface	74
10.3	Reusable Secretion Canister System	75
10.4	Disposable Secretion Canister System	76
11	Maintenance	78
12	Troubleshooting and Service	79

12.1 ASSKEA proVisio^{Basic} 79

12.2 ASSKEA proVisio^{Basic} battery 80

13 Decommissioning, Storage, Disposal 81

14 Patient change in the Homecare sector 82

15 Technical Data 83

16 Order information..... 84

17 Notes on electromagnetic compatibility 85

17.1 Safety instructions..... 85

17.2 Electromagnetic environment in which the ASSKEA proVisio^{Basic} bzw.
 ASSKEA proVisio^{Basic} battery is intended to be operated 86

17.3 Dealing with electromagnetic interaction..... 87

17.4 Overview of all cables and transducers that can be replaced by the user
 87

en

1 Important Information

DAMAGE DURING DELIVERY

Please check shipments immediately upon delivery for completeness and integrity. All defects must be reported within 3 days, otherwise they cannot be accepted.

CUSTOMER FEEDBACK

The products described in these instructions are subject to continuous development and improvement. We therefore welcome any customer feedback, comments and suggestions regarding our products and accompanying documentation that will help us improve our products, services, or documentation.

Please feel free to share your opinion with us!
You can use one of the following options to do so:

- by email to mps@asskea.de
- by telephone on +49 (0) 36201-5797-0

NOTES ON THIS MANUAL

Please read these instructions for use thoroughly before operating the device for the first time. Observe the safety instructions in these instructions for use to avoid hazardous situations.

If you require assistance with commissioning, operation or maintenance, please contact ASSKEA GmbH.

Please report any unexpected operation or incidents to ASSKEA GmbH.

Keep these operating instructions in an easily accessible place.

Please include these operating instructions when passing on the device to third parties.

ILLUSTRATIONS AND SYMBOLS IN THESE INSTRUCTIONS FOR USE

Illustration	Meaning
	Caution – Information on risks and their avoidance.
	Note with useful information and tips
•	List
1. 2.	Perform the action in the order described.

2 Intended Purpose

INTENDED USE

The **ASSKEA proVisio^{Basic}** and **ASSKEA proVisio^{Basic} battery** are intended for temporary removal of suction material in the head and neck region from natural or artificial body openings.

INTENDED PLACE OF USE

- Upper respiratory tract (nose, nasal cavity, throat) and lower respiratory tract (larynx, trachea, bronchial system) for suction of body fluids (blood, secretions) or food components.
- External ear canal to remove cerumen.

en

INDICATIONS

Damage to the respiratory and cough functions with disturbance of the tracheal, bronchial, or oral secretion elimination:

- Tracheostomy patients
- Laryngectomy patients
- Patients with impaired respiratory function

Muscular and/or neurological diseases:

- Patients with swallowing disorders
- Patients with disorders leading to accumulation of blood, secretions and food components from the oral cavity, the pharynx or the bronchial system

Adult patients with a cerumen impaction (cerumen plug) in the external ear canal with present symptoms or if the cerumen interferes with an ENT examination.

INTENDED PATIENT GROUP

All people without restrictions, that fall under the patient selection criteria.

PATIENT SELECTION CRITERIA

- Patients who would benefit from upper or lower airway suction.
- Adult patients who would benefit from removal of cerumen from the external ear canal.

CONTRAINDICATIONS

The decision to use the device must be made from a medical perspective, taking into account the patient's general condition.

Absolute contraindications to nasopharyngeal/oropharyngeal suctioning:

- Epiglottitis or croup

The risks of suctioning body fluids (blood, secretions, food components) from the upper respiratory tract (nose, nasal cavity, throat), lower respiratory tract (larynx, trachea, bronchial system) should be particularly considered in patients with the following symptoms:

Relative contraindications to nasopharyngeal/oropharyngeal suctioning:

- Unexplained hemoptysis or hemoptysis (coughing up blood from the lungs or bronchi) or known coagulopathy
- Laryngospasm
- Bronchospasm
- Basal skull fractures and other causes of spinal fluid leakage from the ear
- Pneumothorax
- Fresh anastomosis of the esophagus or trachea or other tracheobronchial trauma
- Obstructed nasal airway
- Nasal bleeding of unknown cause
- Severe hypoxemia/hypoxia
- Increased intracranial pressure
- Acute hypo- or hypertension
- Recent acute cardiac insult/ event
- High anastomotic surgery (tracheal or oesophageal/ gastric)
- Occluded nasal passages with a history of nosebleeds
- Upper respiratory tract infection (for fear of transfer of infection to lower respiratory tract by the procedure itself)

Relative contraindications to oral and orotracheal suctioning:

- Facial fractures, acute head, face, or neck injury (especially skull fractures or leakage of spinal fluid from the ear)
- Loose teeth or crowns
- (Serious) coagulopathy or hemoptysis
- Laryngeal or oral carcinoma
- Severe bronchospasm
- Laryngospasm / stridor
- Restless / anxious patient
- Fresh anastomosis of the esophagus or trachea or tracheoesophageal fistula
- Gag reflex

Relative contraindications to nasotracheal suctioning:

- Occluded nasal passages
- Nasal bleeding
- Acute head, facial, or neck injury
- Coagulopathy or bleeding disorder
- Laryngospasm
- Irritable airway

- Upper respiratory tract infection
- Tracheal surgery
- Gastric surgery with high anastomosis
- Myocardial infarction
- Bronchospasm
- Recent surgery of the face or nose
- Recent anastomosis of the esophagus or trachea or tracheoesophageal fistula

There is no contraindication to endotracheal suction.

Removal of cerumen is contraindicated in patients with the following symptoms:

- Hyperacusis or tinnitus
- Uncontrollable head movements
- Trauma to the ear canal or eardrum
- Anomalies that prevent a clear view of the application
- Diabetes mellitus
- anticoagulated or immunocompromised state
- prior radiation therapy to the head or neck affecting the ear canal
- If the patient feels pain or discomfort, the treatment must be stopped immediately.

en

RISKS AND SIDE EFFECTS

For the indication

- Damage to the respiratory and cough functions with disturbance of the tracheal, bronchial, or oral secretion elimination
- Muscular and/or neurological diseases

the following risks were identified:

Slow vacuum build-up and insufficient flow can prolong suction time. This can lead to temporarily reduced oxygen saturation or hypoxia/hypoxemia.

To maintain the performance of the device, maintenance procedures must be performed as described in the user's manual.

Contaminated suction devices can introduce pathogens into the respiratory tract. The device does not come into direct contact with the patient.

However, to reduce the general risk of contamination, the device must be cleaned and disinfected as described in the user's manual.

Other non-device-related risks and side effects may occur during suctioning. It is recommended that the user be aware of the general risks associated with airway treatment and seek medical advice promptly in the event of unusual symptoms or complications.

For the indication

- Patients with a cerumen impaction (cerumen plug) in the external ear canal with present symptoms or if the cerumen interferes with an ENT examination.
- the following risks were identified:

Too loud noise can cause tinnitus or hearing loss.

Use the lowest vacuum setting possible to minimize suction noise.

Contaminated suction devices can introduce pathogens into the respiratory tract. The device does not come into direct contact with the patient. **However, to reduce the general risk of contamination, the device must be cleaned and disinfected as described in the user's manual.**

WHO MAY OPERATE THE DEVICE?

- **medical staff, professional user**
- **non-medical staff, e.g. patients or patient's relatives** who have been instructed and trained in handling, operation, reprocessing, and disposal of the device and its accessories in accordance with the instructions for use. The instruction and training is provided by the trained and instructed medical staff/professional user.

Use on the external ear canal is excluded for this user group.

IN WHICH ENVIRONMENT MAY THE DEVICE BE USED?

The device can be operated in:

- in outpatient and inpatient care (professional healthcare facility environment)
- in homecare, especially for aspiration in tracheostomy patients.

The device must not be operated:

- in medical rooms where potential equalization is required (e.g. cardiac surgery).
- in areas at risk of explosion.
- in the MRI environment.
- when used outdoors
- during transport
- in areas with high humidity or damp rooms (e.g. near inhalers or kettles, when bathing or showering).
- in areas subject to excessive heat. (e.g. directly next to heat sources or areas exposed to direct sunlight).

OPERATING MODES

ASSKEA proVisio^{Basic}:	Mains-powered, stationary operation
ASSKEA proVisio^{Basic} battery:	Mains-independent, mobile operation

3 Safety Instructions

WARNING OF PERSONAL INJURY DUE TO IMPROPER USE

- For your own protection and the safety of patients, never use the device for purposes other than those for which it is intended!
- Do not use the device if there is visible damage to the housing.
- Do not open the device if you are not trained or authorised to do so. Only authorised service personnel should carry out service activities such as retrofitting a battery or replacing a battery.
- Modifications to the device are not permitted.

WARNING OF HAZARD DUE TO IMPROPER HANDLING IN HOME USE

Children and pets may unintentionally operate the device or knock it over. This may lead to malfunction or injury caused by falling parts.

- Keep children and pets away from the device and do not leave the device unattended within their reach.

WARNING OF MALFUNCTION DUE TO ENVIRONMENTAL CONDITIONS IN HOME USE

Dust, lint or pet hair may impair the function of the device.

- Keep the device clean and free from dust and lint.

Direct sunlight or heat sources may lead to overheating of the device.

- Do not operate the device in the immediate vicinity of heat sources and protect it from direct sunlight.

High humidity, e.g. caused by inhalers or kettles, may enter the device and lead to malfunction or electrical damage.

- Keep the device away from inhalers or kettles.

WARNING CONCERNING INJURIES DURING TREATMENT

- Select the suction setting in depending on the patient and the medical indication to avoid possible personal injury.
- Perform the applications carefully to avoid injuries in the suction area.
- Never use the device to suction solid objects or inhaled foreign bodies from the throat. This may cause the object to be pushed deeper into the airways, where it may cause total obstruction.
- To ensure patient care in the event of device failure, we recommend having a replacement device available. Call an ambulance if there is a risk of suffocation.

WARNING OF HAZARD DUE TO MISSING DISCONNECTION FROM THE MAINS SUPPLY

The device is connected to the mains supply via the power supply unit. The mains plug serves as the disconnection device from the mains supply. In hazardous situations, the device must be able to be disconnected quickly from the mains supply.

- Position the device so that the mains plug is always easily accessible.
- The mains plug must not be covered, obstructed or difficult to reach, so that disconnection from the mains supply can be carried out at any time.

WARNING AGAINST INJURIES CAUSED BY STRANGULATION

People can strangle themselves with tubes or power cords.

- Ensure that unauthorised persons or bystanders do not remain in the vicinity of the device.
- Remove cables and tubes from the patient's environment when they are no longer needed.

WARNING OF HEALTH RISKS WHEN HANDLING INFECTIOUS OR PATHOGENIC GERMS

Infectious and pathogenic germs in the suctioned material cause health hazards.

- Never operate the device without a bacteria filter.
- Always aspirate using a sterile disposable catheter or a sterile cannula.
- Ensure that the suction tube does not come into direct contact with the suction site.
- Wear disposable gloves to avoid coming into contact with the suctioned material.
- Clean the device and its accessories as described in these instructions for use.

WARNING REGARDING UNAUTHORISED ACCESSORIES AND SPARE PARTS

- Only use the accessories and spare parts listed in these instructions for use. The safety and functionality of the device cannot be guaranteed for other accessories.
- Only use clean accessories that are free of damage.

WARNING – DEVICE DAMAGE DUE TO INGRESS OF LIQUIDS OR SOLIDS

- Keep the device clean and dry.
- Clean the device only as described in this manual.

If liquids or solids enter the interior of the device, it must be inspected by ASSKEA GmbH or an authorized service partner.

CAUTION: DAMAGE TO THE DEVICE DUE TO INCORRECT MAINS CONNECTION

To avoid damage to the device, please observe the following instructions:

- Only use the power supply unit supplied (GTM46402-3713.8).
- When using the power supply unit, make sure that the power supply unit plug is connected to the device first and only then is the power supply unit connected to the mains.
- The power supply unit must be disconnected from the mains connection in exactly the opposite order.

CAUTION: RISK OF DAMAGE TO THE DEVICE DUE TO INCORRECT HANDLING

- Never extract flammable, corrosive or explosive liquids or gases.
- Always place the device upright on a firm, level surface and ensure that it cannot be knocked over or fall down, so that people cannot be hit by the falling device.

CAUTION: DAMAGE TO THE DEVICE DUE TO HEAT BUILD-UP

- Do not cover the power supply unit.
- Keep the suction device, power cord and power supply away from other heat sources.

BEWARE OF RADIO INTERFERENCE.

- Avoid using this device immediately next to other devices or stacked with other devices, as this could result in malfunction. If use in the manner described above is nevertheless necessary, this device and the other devices should be monitored to ensure that they are functioning properly.

The device complies with the requirements of standard EN 60601-1-2. The system is not affected by electromagnetic radiation from other devices that comply with the same standard. Furthermore, the system does not generate electromagnetic radiation that exceeds the limits specified in EN 60601-1-2. Further information on electromagnetic immunity and electromagnetic radiation is available in chapter 17.

en

4 Warranty and Product Liability

WARRANTY

The warranty period for the device is 2 years. It is neither extended nor renewed by warranty work carried out. The warranty period for the battery is 6 months.

Wear parts are excluded from the warranty.

All warranty claims shall expire if:

- the device is opened by unauthorised persons,
- the security seal is removed/damaged,
- repairs are carried out by unauthorised persons,
- changes are made to the device,

as in these cases the basic safety of the device can no longer be guaranteed, and functional restrictions may occur.

Any warranty is excluded for damage caused using non-recommended accessories and spare parts or improper use.

ASSKEA GmbH guarantees service for the device for at least the product life of 5 years.

PRODUCT LIABILITY

Liability for the function of the device is transferred to the operator if:

- the device is used outside its intended purpose,
- the device is not used in accordance with the instructions for use,
- the device is opened by unauthorised persons,
- the safety seal is removed/damaged,
- assembly, adjustments, extensions, modifications or repairs are carried out by unauthorised persons,
- non-approved accessories and spare parts are used.

5 Symbols and Terms

5.1 Symbols on Device, Packaging and Accessories



Observe the symbols on the device, the packaging, and the accessories.

SYMBOLS THAT IDENTIFY THE MEDICAL DEVICE AND ACCESSORIES:

Symbol	Meaning	Symbol	Meaning
	Manufacturer		Date of manufacture with country of manufacture
	Distribution partner		Importer
	Swiss Authorised Representative		The product is a medical device.
	Catalogue number		Serial number
	Batch number		Unique device identification

en

INSTRUCTIONS FOR HANDLING THE MEDICAL DEVICE AND ITS ACCESSORIES:

Symbol	Meaning	Symbol	Meaning
	Follow the instructions in the user manual.		Observe the instructions for use.
	The product is intended for reuse on a single patient.		Protection class: IP22 : The device is protected against finger access and dripping water at an angle of up to 15°.
	The product is intended for single use only.		
	Do not use if the packaging is damaged!		
	Operating time		This is a Class II electrical device. The device is not connected to the protective earth conductor.
	Power supply		
	The application part of the device is classified as type BF . (Body Floating)		This device must not be disposed of with household waste.

INFORMATION ON ENVIRONMENTAL, STORAGE AND TRANSPORT CONDITIONS:

Symbol	Meaning	Symbol	Meaning
	Protect from moisture.		Observe the temperature limits during operation, storage or transport.
	Observe the air pressure limit during operation, storage or transport.		Observe the air humidity limit during operation, storage or transport.

5.2 Glossary

A

Applied part type BF Applied parts are classified according to their area of application.
Type BF application parts must be operated isolated from earth and are not suitable for direct application to the heart.

C

Contamination Contamination means that bacteria or viruses from the suction material have come into contact with the device.

O

Over-suction Over-suction means that the suctioned material is sucked into the interior of the device.

S

Suction material Suction material is a generic term for body fluids (blood, secretions), food residues, or cerumen. These can be easily suctioned off using the device described here.

6 Product Description

6.1 Device and Accessories - Overview

DEVICE VARIANTS

ASSKEA proVisio^{Basic}: without battery, for stationary operation
 ASSKEA proVisio^{Basic} battery: with battery, for mobile operation

DEVICE COMPONENTS AND ACCESSORIES

Both device variants can be operated with a **reusable secretion canister system** or a **disposable secretion container system**.

en

DEVICE WITH REUSABLE SECRETION CANISTER SYSTEM


Fig. 1

- 1 Carrying handle
- 2 Vacuum gauge with vacuum control knob
- 3 Display and control panel
- 4 Power supply socket
- 5 Power supply unit
- 6 visioBag (Reusable secretion canister with lid and tube mount)
- 7 Suction tube with fingertip
- 8 Rinsing bottle

DEVICE WITH DISPOSABLE SECRETION CANISTER SYSTEM


Fig. 2

- 1 Carrying handle
- 2 Vacuum gauge with vacuum control knob
- 3 Display and control panel
- 4 Power supply socket
- 5 Power supply unit
- 6 Disposable secretion canister and Holder for external canister "Bag"
- 7 Suction tube with fingertip
- 8 Rinsing bottle

REUSABLE SECRETION CANISTER SYSTEM

Fig. 3: Reusable secretion canister


- 1 Lid for reusable secretion canister
- 2 reusable secretion canister

Fig. 4: Mechanical overflow protection

- 1 Overflow protection with ball

Fig. 5: Bacteria filter of the reusable secretion canister system

- 1 Bacteria filter of the reusable secretion canister system


DISPOSABLE SECRETION CANISTER SYSTEM

Fig. 6: Disposable secretion canister


- 1 Disposable liner "OneWay"
- 2 External canister "Bag"
- 3 T-piece
- 4 Connecting tube

Fig. 7: Connecting ports of disposable liner "OneWay"

- 1 Vacuum connection
- 2 Patient connection

Fig. 8: Adapter for the disposable secretion canister system

- 1 Adapter for the disposable secretion canister system

Fig. 9: Holder for external canister "Bag"

- 1 Holder for external canister "Bag"

6.2 Scope of Delivery

STANDARD PACKAGE CONTENTS

- **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery**, including
 - Instructions for use
 - Power supply unit (GTM46402-3713.8)
- visioBag (Item no. 100281)
- Suction tube with fingertip (Item no. 100713-3)
- Rinsing bottle (Item no. 100282-1)
- Additional components may be included depending on the order.

en

6.3 Applied Parts

In order to operate the device for its intended purpose, a type BF application part is also required. The following application parts can be connected to the suction tube:

- Suction catheter with conical plug connection
- Suction cannula with Luer Lock for use in the external auditory canal

Select the application part depending on the patient, the medical indication and the viscosity of the suction material.

When using this suction device at home, the attending physician should determine the appropriate size of the catheter. Please consult your physician to ensure that the catheter used is correctly matched to your individual needs.

If the catheter or cannula diameter is too small, suction may not be performed adequately.

Effective suctioning of newborns and paediatric patients with impaired respiratory and cough function with impaired tracheal, bronchial or oral secretion clearance and patients with muscular and/or neurological disorders is possible with a suction catheter size of CH6 or larger; in older patients, this is possible with CH8 or larger.

For use in the external auditory canal to remove cerumen, use 17G cannulas or cannulas with a larger diameter.

6.4 Connection between Device and Accessories

The diagram below shows how the device is connected to the accessories and patient.



Fig. 10:

Assembling the device with the secretion canister system and the suction tube is described in Chapter 7.


You can connect the suction tube to the various application parts as follows.

Suction tube with fingertip → Suction catheter with conical connection

Connection	Type	Diameter / Standard
Suction tube → Fingertip	Conical plug connection	6–10 mm / non-standard
Fingertip → Suction catheter	Conical plug connection	6–10 mm / ISO 5356-1

Suction tube with fingertip → adapter → suction cannula with Luer lock for use in the external auditory canal

Connection	Type	Diameter / Standard
Suction tube → Fingertip	Conical plug connection	6–10 mm / non-standard
Fingertip → Adapter	Conical plug connection	6–10 mm / non-standard
Adapter → Suction cannula with Luer-Lock	Luer lock connection	Ø 4 mm / ISO 80369-7

6.5 Battery

CHARGE THE BATTERY

The battery is not automatically kept in an operational state. Charge the device when the battery LED status indicator lights up orange.

The **ASSKEA proVisio^{Basic} battery** is equipped with a lithium-ion battery. The battery is protected by safety measures against deep discharge and overheating during charging.

If, due to improper ambient conditions, the battery temperature exceeds the permissible level during charging, the charging process will be temporarily interrupted to allow the battery to cool down. This measure ensures safety and helps to protect the battery.

The typical service life of the battery is approximately 500 charging cycles, after which around 70% of the original capacity remains.

PLEASE OBSERVE THE FOLLOWING INSTRUCTIONS:

- Fully charge the device before first use.
- Repeat the charging process after the first few uses.
- Avoid frequent short charging cycles.

en

6.6 Reusable Secretion Canister System

The components of the reusable secretion canister system are intended for reuse on a single patient.

OVERFLOW PROTECTION

The lid of the reusable secretion canister system contains overflow protection (float shut off). It consists of a float cage and a float ball integrated into the lid. When the canister is full, the ball rises within the cage and closes, in the upper section, the passage from the lid to the device. This interrupts the suction process and prevents the device from over-suctioning.

BACTERIA FILTER

The bacteria filter between the device and the lid of the reusable secretion canister is made of a hydrophilic (self-sealing) material and is effective against bacteria and viruses. It effectively protects the device from over-suction and contamination.

Please observe the replacement intervals in Chapter 11.

6.7 Disposable Secretion Canister System

The components of the disposable secretion canister system are intended for reuse on a single patient.

OVERFLOW PROTECTION / BACTERIA FILTER

The disposable liner "OneWay" has an integrated bacteria filter, carbon filter and gelling agent. The self-sealing bacteria filter integrated into the disposable liner "OneWay" acts against microorganisms. In the event of an application error, this integrated filter prevents over-suction. When liquid reaches this filter, it blocks and further suction is no longer possible. The suction process is interrupted. The disposable liner "OneWay" must be replaced. The activated carbon filter in the disposable liner "OneWay" reduces the spread of odours.

GELLING AGENT

Disposable liner "OneWay" filled with suction material can be transported and disposed of in a leak-proof manner with the aid of the gelling agent. Regardless of the suction intervals, the suction material is thickened after an average gelling time of 2-5 minutes (depending on the nature of the secretion).

The **adapter for the disposable secretion canister system** is marked "no filter" (Fig. 8) and may only be used with the disposable secretion canister system, as the adapter itself does not have a bacteria filter.

The bacteria filter is in the lid of the disposable liner "OneWay". This effectively protects the vacuum generator from over-suction.

Please observe the replacement intervals in Chapter 11.

7 Start-up

7.1 Starting up the device

CHECK POWER SUPPLY UNIT AND MAINS CABLE

Check the power supply unit, including mains cable for possible damage or defects and replace it immediately if visible damage is present.

Allow the device to acclimatise at room temperature (approx. 20 °C) for about 2 hours before switching it on to prevent damage.

en

ASSKEA proVisio^{Basic}

1. Remove all packaging from the device and accessories.
2. Place the device on a firm, level surface.
3. Connect the power supply plug to the suction unit.
4. Connect the power supply to the mains.
5. Connect the reusable secretion container system (Section 7.2) or the disposable secretion container system (Section 7.3) with the device.

ASSKEA proVisio^{Basic} battery

1. Remove all packaging from the device and accessories.
2. Place the device on a firm, level surface.
3. Connect the power supply plug to the suction unit.
4. Connect the power supply to the mains.
5. Fully charge the battery.
6. Connect the reusable secretion container system (Section 7.2) or the disposable secretion container system (Section 7.3) with the device.

7.2 Connecting the Reusable Secretion Canister System

CHECK BACTERIA FILTER

Ensure that the bacteria filter is inserted in the device. Always have an additional bacteria filter ready, as this is essential for safe operation.

CHECK REUSABLE SECRETION CANISTER

Before each use, check that the reusable secretion canister is completely clean to prevent foaming.

1. Fold the movable part of the secretion canister lid (visioClick®) slightly upwards.


Fig. 11

2. Place the secretion canister lid evenly on the secretion canister.


Fig. 12

3. Fold the movable part of the lid (visioClick®) downwards so that the side hooks engage and the secretion canister is securely closed.


Fig. 13

4. Insert the secretion canister into the designated canister holder on the device by guiding it over the bar at the bottom


Fig. 14

5. Then click the upper part of the secretion container into place.


Fig. 15

6. Connect the suction tube to the connection in the middle of the secretion canister lid. Ensure that the tube is securely attached.


Fig. 16

7.3 Connection of the Disposable Secretion Canister System

Malfunction due to collapsing disposable liner "OneWay"

Air can flow into the external canister "Bag" through a leak in the external canister or in the lid of the disposable liner "Oneway". This can cause the disposable liner "Oneway" to collapse.

- Ensure that the lid of the disposable liner "Oneway" is firmly connected to the external canister "Bag".
- Ensure that all connections are secure and correctly connected.
- Ensure that the external canister "Bag" is not damaged and that the T-piece is secure.

Always have an additional disposable liner "OneWay" ready, as this is essential for safe operation!

en

1. Ensure that the adapter for the disposable secretion canister system is correctly inserted. After insertion, rotate the adapter handle a quarter turn to the left until the lever is horizontal.


Fig. 17


Fig. 18

2. Position the holder for the external canister "Bag" at the top, under the notch of the canister mount.
3. Then slide it at the bottom over the bar fully into the device.


Fig. 19


Fig. 20

4. Remove the disposable liner "OneWay" from the packaging and fully unfold it.
5. Insert it into the external canister "Bag" and press the lid firmly down along the rim to ensure a leak-tight seal.


Fig. 21


Fig. 22

6. Insert the lower end of the T-piece of the external canister "Bag" into the opening in the holder. Make sure that that the connection is correct.

7. Attach the pre-assembled connecting tube of the disposable liner "OneWay" to the upper end of the T-piece.


Fig. 23


Fig. 24

8. Connect the suction tube to the patient connection of the disposable liner "OneWay" ②. Ensure the tubing is firmly seated.



Fig. 25

- ① Vacuum connection
- ② Patient connection

Fig. 26

8 Operation

8.1 Control and Display Elements

VACUUM GAUGE WITH VACUUM CONTROL KNOB

Fig. 27:

Vacuum setting:

-100 mbar to -850 mbar

-10 kPa to -85 kPa

en

DISPLAY AND CONTROL PANEL


Fig. 28: ASSKEA proVisio^{Basic} battery left, ASSKEA proVisio^{Basic} right

- 1** Power Button

 -  OFF Device is switched OFF
 -  GREEN Device is switched ON
- 2** Power Supply Status LED

 -  OFF No power supply connected
 -  GREEN Power supply is connected
 -  RED Wrong power supply connected
- 3** Battery discharge-status LED

 -  OFF Battery charge level more than 10 %
 -  ORANGE Battery charge level below 10 % or Temperature Error.
- 4** Battery charge-status LED

 -  OFF No power supply connected
 -  ORANGE Power supply is connected and battery is charging
 -  Blinking ORANGE Power supply connected and Battery ERROR!

8.2 Operation and Function

When switching on the device, the pump starts and generates a vacuum within the canister system. This vacuum enables the suction of fluids or material via a sterile suction catheter or a sterile suction cannula for use in the external auditory canal.

The vacuum can be adjusted using the vacuum control knob. Adjustment is continuous up to a maximum vacuum of - 85 kPa or - 850 mbar ($\pm 10\%$). The vacuum is technically limited to -90 kPa. The maximum suction capacity is 24 L/min ($\pm 10\%$).

The device automatically switches off after 5 minutes of operating time, regardless of whether suction is actively performed or the suction tube is open to air.

An overheating protection device also prevents the plug-in battery from overheating by interrupting the charging process if the plug-in battery becomes too hot (e.g. due to unfavourable environmental conditions).

If the device is used in mains operation, it switches off automatically as soon as the power plug is disconnected. There is **no automatic switchover to battery operation**. Similarly, the device switches off automatically in battery operation when the power supply unit is plugged in to charge the battery.

For immediate further use, the device can or must be **switched on again immediately** after changing the power plug.


8.3 Functional Check

1. Switch on the device.
2. Adjust the vacuum using the vacuum control knob.
3. Cover the opening of the suction tube to check whether a vacuum is present.
 - ✓ Vacuum present:
Suction can be performed.
 - ✗ No vacuum present:
Eliminate the cause of the fault according to Chapter 12 or contact the customer service.

8.4 Performing Suction

1. Switch on the device.
2. Adjust the vacuum using the vacuum control knob.
3. Perform the suction procedure.
The suction power can be fine-adjusted using the air regulation on the fingertip control.
4. Dispose of the applied part afterwards.
5. Rinse the suction tube with water after use to prevent it from drying out. Ensure that the secretion container has sufficient capacity!
6. Press the Power button to switch off the device.

en

9 Follow-up

WARNING OF HEALTH HAZARDS WHEN HANDLING INFECTIOUS OR PATHOGENIC GERMS

Infectious and pathogenic germs in the suctioned material cause health hazards.

Wear disposable gloves to avoid coming into contact with the suctioned material.

9.1 Reusable secretion canister system is full and must be emptied

1. Switch off the device after suctioning.
2. When operating from the mains, first unplug the mains plug from the wall outlet, then disconnect the device plug from the suction unit.
3. Remove the suction tube from the application part.
4. Remove the suction tube from the secretion canister.
5. Press the release mechanism on the canister lid to separate the canister from the device.


Fig. 29


Fig. 30

- Remove the reusable secretion canister from the device.


Fig. 31

- Fold the movable part of the lid (visioClick®) upward to release the side hooks, then remove the lid.

Fig. 32

- Empty and clean the reusable secretion canister and lid or dispose of the secretion canister system when changing patients.

9.2 Disposable secretion canister system is full and must be replaced

- Switch off the device after suctioning.
- When operating from the mains, first unplug the mains plug from the wall outlet, then disconnect the device plug from the suction unit.
- Remove the suction tube from the application part.
- Remove the suction tube.
- Loosen the pre-assembled connecting tube of the disposable liner "OneWay" at the upper end of the T-piece of the external canister "Bag".


Fig. 33


Fig. 34

6. Remove the external canister "Bag" with the suction bag "OneWay" from the device.
7. Remove the disposable liner "OneWay" from the reusable external canister "Bag".


Fig. 35


Fig. 36

8. Dispose of the disposable liner "OneWay".
9. Insert a new disposable liner "OneWay". See chapter 7.3 (steps 4–8).
10. Dispose of reusable components of the disposable secretion canister system when changing patients.

en

10 Cleaning and Disinfection

10.1 WARNINGS

WARNING OF HEALTH HAZARDS WHEN HANDLING INFECTIOUS OR PATHOGENIC GERMS

Infectious and pathogenic germs in the suctioned material can cause damage to health.

- Wear at least suitable disposable gloves to avoid coming into contact with the suction material.
Depending on the patient's illness, additional protective measures may be necessary. Observe the operator's specifications.
- Carefully follow all the steps listed below to achieve the required reduction in germs.

WARNING OF POSSIBLE PHYSICAL INJURY DUE TO ELECTRIC SHOCK

- Switch off the device before cleaning.
- Disconnect the power supply unit from the mains and then disconnect the power supply unit from the device.

CAUTION: HEALTH HAZARDS WHEN HANDLING DISINFECTANTS

- Follow the manufacturer's instructions for the disinfectant.

BEWARE OF DAMAGE TO THE DEVICE DUE TO INCORRECT CLEANING AGENTS

Cleaning agents and disinfectants can damage the surface of the device.

- Do not use aldehydes or acetone for cleaning and disinfection.
- Frequent cleaning and disinfection may cause slight discolouration of the plastic parts of the housing. However, this does not affect the functionality of the device.

10.2 Device Surface

Clean and disinfect the surfaces of the device regularly, at least once a week and between patients.

Limitations on reprocessing	Do not use the device if it is visibly damaged or if suctioned material has entered the device. There are no further restrictions.		
Preparation before cleaning	Follow the instructions in Chapter 9. The device should now be disconnected from the mains.		
Cleaning	Remove visible dirt with a damp, but not wet, lint-free cloth. Moisten the cloth with clean water (< 40°C). Wipe the entire surface of the device with light pressure.		
Disinfection	Now disinfect the previously cleaned surfaces. Follow the manufacturer's instructions for the disinfectant. Spray the entire surface of the device. Wipe the entire surface. The device can be disinfected with the disinfectant group "alcohol". ASSKEA GmbH recommends the following disinfectants:		
	Surface disinfectants	Active ingredient group	Type
	Incidin® Liquid (Ecolab)	Alcohol	Ready-to-use solution
	Bacillol® AF (Bode Chemie)	Alcohol	Ready-to-use solution
Drying	Allow the disinfected surface to dry completely.		
Maintenance, inspection and testing	The device should now be free of visible contamination. You can now prepare the device again according to Chapter 7 for the next application.		
Storage	If you are not going to use the device for a longer period of time, store it according to the instructions in Chapter 13.		

10.3 Reusable Secretion Canister System

Clean and disinfect the surfaces of the reusable components of the reusable secretion canister system regularly.

- Reusable secretion canister: as needed, at least weekly
- Lid of reusable secretion canister: as needed, at least weekly
- Suction tube: as needed

The reusable secretion canister, the lid of the reusable secretion canisters and the suction tube must be replaced when changing patients!

en

<p>Limitations on reprocessing</p>	<p>Replace the reusable secretion canister and the lid of the reusable secretion canisters:</p> <ul style="list-style-type: none"> - if damaged, - if cleaning is no longer possible due to heavy soiling, or - when changing patients. <p>Replace the suction tube:</p> <ul style="list-style-type: none"> - if it can no longer be cleaned due to heavy soiling, - at least once a week, or - when changing patients.
<p>Preparation before cleaning</p>	<p>Remove the reusable secretion canister system and separate the components of the system as described in Chapter 9.</p>
<p>Cleaning</p>	<ol style="list-style-type: none"> 1. Empty the reusable secretion canister and rinse it thoroughly with clean water (< 40°C). Use a lint-free cloth and apply light pressure to remove all visible contamination. 2. Remove the overflow protection, including the ball, from the secretion canister lid. <div data-bbox="587 1507 916 1774" data-label="Image"> </div> <ol style="list-style-type: none"> 3. Rinse all parts thoroughly with clean water (< 40°C) as described above for the reusable secretion canister. 4. Rinse the suction tube thoroughly with clean water (< 40°C).

	The reusable secretion canister, the suction tube and the lid for reusable secretion canisters (without overflow protection) are dishwasher-safe and can also be boiled after rinsing.		
Disinfection	You can disinfect the previously cleaned surfaces of the accessories. Follow the manufacturer's instructions for the disinfectant.		
	The accessories can be disinfected with the disinfectant group "alcohol".		
	ASSKEA GmbH recommends the following disinfectants:		
	Surface disinfectants	Active ingredient group	Type
	Incidin® Liquid (Ecolab)	Alcohol	Ready-to-use solution
	Bacillol® AF (Bode Chemie)	Alcohol	Ready-to-use solution
Drying	Allow the cleaned and disinfected surfaces to dry completely.		
Maintenance, inspection and testing	The accessories should now be free of visible dirt. You can prepare the device for the next use again in accordance with Chapter 7.		
Storage	Store the parts in a clean and dry place until the next use.		

10.4 Disposable Secretion Canister System

Clean and disinfect the surface of the reusable components of the disposable secretion collection system regularly.

- Holder for the external canister "Bag": as needed, at least weekly
- External canister "Bag": as needed, at least weekly
- Disposable liner "OneWay" as needed
- Suction tube: as needed

The external canister "Bag", the holder for the external canister "Bag", the disposable liner "OneWay" and the suction tube must be replaced when changing patients!

Limitations on reprocessing	<p>Replace the external canister "Bag" and the holder for the external canister "Bag":</p> <ul style="list-style-type: none"> - if damaged, - if cleaning is no longer possible due to heavy soiling, or - when changing patients.
------------------------------------	---

	<p>Replace the disposable liner "OneWay":</p> <ul style="list-style-type: none"> - when it is full, - at least once a week, or - when changing patients. <p>Replace the suction tube:</p> <ul style="list-style-type: none"> - if it can no longer be cleaned due to heavy soiling, - at least once a week, or - when changing patients. 									
Preparation before cleaning	<p>Remove the disposable secretion canister system and separate the components of the system as described in Chapter 9.</p> <p>Dispose of the disposable liner "OneWay".</p>									
Cleaning	<ol style="list-style-type: none"> 1. Clean visible soiling on the holder external canister "Bag" and the external canister „Bag" with a damp, lint-free cloth or under running clear water (< 40°C) using a lint-free cloth. 2. Rinse the suction tube thoroughly with clean water (< 40°C). <p>The external canister "Bag" and the suction tube is dishwasher-safe and can also be boiled after washing.</p>									
Disinfection	<p>You can disinfect the previously cleaned surfaces of the accessories.</p> <p>Follow the manufacturer's instructions for the disinfectant. These can be disinfected with the disinfectant group "alcohol". ASSKEA GmbH recommends the following disinfectants:</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Surface disinfectants</th> <th style="text-align: left;">Active ingredient group</th> <th style="text-align: left;">Type</th> </tr> </thead> <tbody> <tr> <td>Incidin® Liquid (Ecolab)</td> <td>Alcohol</td> <td>Ready-to-use solution</td> </tr> <tr> <td>Bacillol® AF (Bode Chemie)</td> <td>Alcohol</td> <td>ready-to-use solution</td> </tr> </tbody> </table>	Surface disinfectants	Active ingredient group	Type	Incidin® Liquid (Ecolab)	Alcohol	Ready-to-use solution	Bacillol® AF (Bode Chemie)	Alcohol	ready-to-use solution
Surface disinfectants	Active ingredient group	Type								
Incidin® Liquid (Ecolab)	Alcohol	Ready-to-use solution								
Bacillol® AF (Bode Chemie)	Alcohol	ready-to-use solution								
Drying	<p>Allow the cleaned and disinfected surfaces to dry completely.</p>									
Maintenance, inspection and testing	<p>The accessories should now be free of visible dirt. You can prepare the device for the next use again in accordance with Chapter 7.</p>									
Storage	<p>Store the parts in a clean and dry place until the next use.</p>									

11 Maintenance

DO NOT OPEN THE DEVICE

Do not open the device if you are not trained or authorised to do so. Only authorised service personnel may carry out service activities.

Keep the device in a generally clean condition.

REPLACEMENT INTERVALS

Bacterial filter of the reusable secretion canister system

- if it has come into contact with suctioned material
- when changing patients

Adapter of the disposable secretion canister system

- if damaged
- when changing patients

Reusable secretion canister, lid of the reusable secretion canister, External canister "Bag" and holder for external canister "Bag":

- if damaged,
- if cleaning is no longer possible due to heavy soiling, or
- when changing patients.

Disposable liner "OneWay":

- when it is full,
- at least once a week, or
- when changing patients.

Suction tube:

- if it can no longer be cleaned due to heavy soiling,
- at least once a week, or
- when changing patients.

Applied part (single-use item)

- dispose after use

REPLACING THE BACTERIAL FILTER AND ADAPTER

1. Remove the bacterial filter or adapter by turning the lever—starting from the horizontal position—a quarter-turn counterclockwise (CCW) until the lever points vertically downward.

Fig. 37


Fig. 38

2. Pull the bacterial filter or adapter out of the socket in this position.
3. Insert the new bacterial filter or adapter in reverse order.


Fig. 39

12 Troubleshooting and Service

en

WARNING - DEVICE DAMAGE DUE TO INGRESS OF LIQUIDS OR SOLIDS

- Keep the device clean and dry.
- Clean the device only as described in this manual.

If liquids or solids enter the interior of the device, it must be inspected by ASSKEA GmbH or an authorized service partner.

Contact ASSKEA GmbH or your service partner if a fault cannot be rectified using the measures described here.

12.1 ASSKEA proVisio^{Basic}

Display	Possible Cause	Solution
DEVICE DOES NOT START		
	Power supply unit or power cable not connected.	Connect the power supply unit and power cable.
	Incorrect power supply unit used.	Connect the correct power supply unit.
INSUFFICIENT SUCTION PERFORMANCE		
	Leakage at tube connections or canister.	Check correct connections with accessories.
	Bacterial filter contaminated or moist.	Replace bacterial filter.
NO SUCTION PERFORMANCE		
	Overflow protection in the canister activated.	Empty reusable secretion canister. Replace disposable liner "OneWay".
	Bacterial filter blocked.	Replace bacterial filter.
	Suction tube clogged.	Rinse or replace suction tube.
	Bacterial filter or adapter not inserted.	Insert bacterial filter or adapter.

Display	Possible Cause	Solution
DEVICE SHUTS DOWN		
	The device turns off automatically after 5 minutes.	Switch the device on again.
	The engine has overheated and switched itself off	Switch off the appliance, leave it to cool down, then switch it back on.

12.2 ASSKEA proVisio^{Basic} battery

Display	Possible Cause	Solution
DEVICE DOES NOT START		
	Power supply unit or power cable not connected.	Connect the power supply unit and power cable.
	Incorrect power supply unit used and battery is empty.	Connect the correct power supply unit.
	Battery temperature too high during discharge.	Allow the battery to cool down and restart.
DEVICE DOES NOT CHARGE.		
	Battery temperature too high during charging.	No action required. Charging starts automatically after cooling down.
	Battery defective.	Contact service!
	Input voltage out of allowable range.	Connect the correct power supply, check the power source.
	Input voltage out of allowable range.	Connect the correct power supply, check the power source.
	Battery temperature too high during discharge.	No action required. Charging starts automatically after cooling down.
INSUFFICIENT SUCTION PERFORMANCE		
 Or 	Leakage at tube connections or canister.	Check correct connections with accessories.
	Bacterial filter contaminated or moist.	Replace bacterial filter.

Display	Possible Cause	Solution
NO SUCTION PERFORMANCE		
  Or	     	Overflow protection in the canister blocked. Empty reusable secretion canister. Replace disposable liner "OneWay".
	Bacterial filter blocked.	Replace bacterial filter.
	Suction tube clogged.	Rinse or replace suction tube.
	Bacterial filter or adapter not inserted.	Insert bacterial filter or adapter.
DEVICE SHUTS DOWN.		
  Or	     	Battery temperature too high during discharge. Allow the device to cool down and switch it on again.
  Or	     	The engine has overheated and switched itself off. Switch off the appliance, leave it to cool down, then switch it back on.
 	     	The device turns off automatically after 5 minutes. Switch the device on again. Battery depleted. Connect the power supply to charge the battery.

en

13 Decommissioning, Storage, Disposal

DECOMMISSIONING

1. After suction is completed, switch the device off.
2. First unplug the mains plug from the wall outlet, then disconnect the device plug.
3. Remove the accessories as described in Chapter 9.
4. Dispose of single-use products.
5. Clean reusable accessories according to Chapter 10.
6. Clean the exterior surfaces of the device as described in Chapter 10.
7. Store the device in the shipping carton until the next use.

STORAGE

Store the device only under the following conditions:

- Ambient temperature: **-25 °C to +70 °C**
- Relative humidity: **up to 90%, non-condensing**
- Air pressure: **700 hPa to 1060 hPa**

Fully charge the battery before storage.

If the device is not used for approximately 10 months, the battery must be recharged.

DISPOSAL

Packaging

Recycle the packaging.

Accessories

- Do not reuse single-use products. Dispose of them properly.
- Reusable products: clean and disinfect, then recycle.

Device

- Clean and disinfect the device.
- Do not dispose of the device or the rechargeable battery with household waste.
- Return both to a recycling centre or return them to ASSKEA GmbH or an authorised service partner.

14 Patient change in the Homecare sector

When changing patients in the home care sector, professional hygienic reprocessing by ASSKEA GmbH or a service partner authorised by ASSKEA GmbH is required.

For this purpose, hand over the device to ASSKEA GmbH or to specialist personnel authorised by ASSKEA GmbH.

TRANSPORT/RETURN

ASSKEA GmbH offers its partners and customers fast and professional reprocessing as well as the performance of necessary tests.

- Clean the device before shipping in accordance with the instructions in Chapter 10.
- Enclose the power supply unit and power cable with the device.
- Dispose of disposable products.
- Dispose of accessories that are intended for use on a single patient only.
- Register the return shipment.
The goods return form can be found on our website in the Service section:
Website: www.asskea.de.

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Model Name / REF	ASSKEA proVisio ^{Basic} / REF 100738-1 ASSKEA proVisio ^{Basic} battery / REF 100738-2
Maximum suction flow rate*	24 L/min ± 10% (high flow)
* Effective flow rate measured at the device connection port. Values may vary with altitude, ambient air pressure, and air temperature.	
Adjustable vacuum range	-10 kPa to -85 kPa
Vacuum tolerance (relative to the set value)	±10 %
Maximum vacuum	-90 kPa
Conversion factors	1 kPa ≈ 10 mbar / 1 kPa ≈ 7,5 mmHg
Operating mode (duty cycle)	3x (5 min ON / 5 min OFF) 120 min COOLDOWN
Protection class <i>against electric shock (EN 60601-1)</i> <i>Ingress protection (EN 60529)</i>	Class II equipment IP 22
Dimensions (H × W × D)	213 mm x 210 mm x 124 mm
Device weight	ASSKEA proVisio ^{Basic} : 1,48 kg ASSKEA proVisio ^{Basic} battery: 1,78 kg
Maximum sound emission	67dB(A)
Environmental conditions	
<i>Transport</i>	Temperature: -25 °C to +70 °C relative humidity: up to 90 %, non-condensing Air pressure: 700 hPa to 1060 hPa
<i>Operation</i>	Temperature: 5 °C bis +40 °C relative humidity: 15 % to 90 %, non-condensing Air pressure: 700 hPa to 1060 hPa

en

ELECTRICAL SPECIFICATIONS

Mains power supply (AC adapter)	In: AC 100 – 240 V ~ 50-60 Hz / 1,0 A Out: DC 13,8 V / 2,68 A / 37 W
Permissible input voltage (device)	13,8 V
Maximum load current	2,68 A
Power consumption	70 W

BATTERY SPECIFICATIONS FOR ASSKEA PROVISIO^{BASIC} BATTERY

Rechargeable battery	min. 10,8 V; lithium-ion
Charging time from empty	min. 2 h
Battery pack energy	<80 Wh
Operating time on battery	approx. 45 minutes, depending on motor run time

ACCESSORY SPECIFICATIONS

Reusable canister	Identification	Reusable secretion canister (1 l) Lid Reusable secretion canister	Item no. 100741, with Item no. 100742
	Container volume	1l	
Disposable canister	Identification	External canister "Bag" Disposable liner "OneWay"	Item no. 100000, with Item no. 100001
	Container volume	1l	
Suction tube		Suction tube CH 25, 130 cm, including fingertip Material: PVC (DEHP-free/latex-free)	
Application parts		Sterile suction catheter or sterile suction cannula, Application part type BF	

16 Order information

Not all items listed in this section are covered by EU MDR 2017/745.

Medical devices or accessories for medical devices are marked with the symbol **MD**.

Item number	Description	PU
100741	Reusable secretion canister (1 l) proVisio®	1
100742	Lid for reusable secretion canister proVisio®	1
100743-1	Reusable secretion canister system proVisio® + suction tube	1
100745-1	Universal holder proVisio®	1
100732-1	5-foot roll stand	1
100748	Device holder for roller stand proVisio®	1
100746	Bag proVisio®	1
100282	Rinsing bottle (250 ml)	1
100747	Tube mount proVisio®	5
100753	Car connection cable	1
100713-3	Suction tube CH 25, 130cm (not sterile)	10
100766	Holder rinsing canister (1 l) proVisio®	1
100767	Set rinsing canister (1 l), lid, holder proVisio®	1
100769	Lid for rinsing canister (1 l) proVisio®	1
100770	Holder external canister „Bag“ (1 l) for proVisio®	1
100772-1	Disposable secretion canister system proVisio®	1
100000	External canister „Bag“ (1 l)	1
100002	Disposable liner „OneWay“ (1 l)	60
400516	Bacterial filter proVisio®	10
400510	Adapter for disposable secretion canister system (1 l)	1

17 Notes on electromagnetic compatibility

17.1 Safety instructions

WARNING: Portable and mobile RF communications equipment (including peripheral devices such as antenna cables and external antennas) may affect medical electrical equipment and should therefore not be used closer than 30 cm to any part of the device, including the cables. Otherwise, the performance of the device may be impaired.

WARNING: The use of accessories and spare parts, converters and cables for the device that are not specified or provided by ASSKEA GmbH may lead to increased electromagnetic emissions or reduced electromagnetic immunity of the device, combined with improper operation.

Any warranty is excluded for damage caused by the use of non-recommended accessories and spare parts, transducers and cables or improper use.

Only use original ASSKEA accessories and spare parts!

WARNING: The use of the specified or provided accessories and spare parts, transducers and cables with devices other than the device may lead to increased electromagnetic emissions or reduced electromagnetic immunity.

Any warranty is excluded for damage caused by the use of the specified or provided accessories and spare parts, converters and cables with other devices or improper use. Only use the accessories and spare parts, converters and cables with the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery!**

WARNING: The use of the device immediately next to or stacked with other devices should be avoided, as this could lead to improper operation. If it is necessary to use the device immediately next to or stacked with other devices, the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** and the other devices should be monitored to verify that they are operating as intended in this configuration!

WARNING: The **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** may interfere electromagnetically with other devices, examinations and treatments. For this reason, special attention should always be paid to other devices and parallel examinations or treatments in order to detect any interference as quickly as possible.

The **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** meets the requirements of EN 60601-1-2 "Electromagnetic compatibility – Medical electrical equipment" without deviations or limitations. Electromagnetic interference and interaction are thus reduced to a minimum. Follow the instructions and guidelines provided to maintain the basic safety and essential functions of the device throughout its expected operating life.

17.2 Electromagnetic environment in which the ASSKEA proVisio^{Basic} bzw. ASSKEA proVisio^{Basic} battery is intended to be operated

The **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** is intended for use in an electromagnetic environment as specified below, in which RF interference is controlled. The customer or user of the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** should ensure that it is used in such an environment.

Intended environments of use include professional healthcare facilities and home care. Special environments such as those near RF surgery or MRI, or environments in which the intensity of EMC interference is high, are excluded.

Emission limits		
Conducted and radiated RF emissions	CISPR 11	
Enclosure		
Phenomenon	Test method	Immunity test level
Static electricity discharge	IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge
Radiated RF disturbances	IEC 61000-4-3	28 V/m 80 MHz to 6 GHz 80% AM at 1 kHz
Radiated RF disturbances	according to frequencies and test levels in EN 60601-1-2, Table 9	see EN 60601-1-2, Table 9
Magnetic field at the supply frequency (50 Hz)	IEC 61000-4-8	30 A/m
Magnetic fields in the immediate vicinity	IEC 61000-4-39	30 kHz, 8 A/m, CW; 134.2 kHz, 65 A/m, PM: 2.1 kHz; 13.56 MHz, 7.5 A/m, PM: 50 kHz
Mains connection AC		
Phenomenon	Test method	Immunity test level
Fast transient electrical disturbances / bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Surge voltages / Surges	IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-earth
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

Phenomenon	Test method	Immunity test level
Voltage dips	IEC 61000-4-11	0% U_T for half a period At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T for one period 70% U_T for 25 periods single-phase at 0°
Voltage interruptions	IEC 61000-4-11	0% U_T for 250 periods

Note: U_T is the mains AC voltage before the test levels are applied.

en

17.3 Dealing with electromagnetic interaction

Although the electromagnetic interference and interaction of **the ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** has been reduced to a minimum, electromagnetic interference between the device and other devices cannot be ruled out. Therefore, always comply with the specified requirements and instructions regarding the permissible electromagnetic environment and monitor the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** to ensure proper functioning and prevent undesirable events for patients and users. If necessary, choose a different location for the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** if the permissible electromagnetic environment cannot be ensured or if you have noticed that the device or other nearby devices are not functioning as intended.

Since the propagation of electromagnetic fields is influenced by absorption and reflection from buildings, objects and people, and the electromagnetic field strength from fixed transmitters cannot be precisely predicted theoretically, a study of the electromagnetic phenomena at the location where the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery** is to be operated in order to determine the electromagnetic environment at the location. If unusual performance characteristics are observed, additional measures may be necessary, such as changing the orientation or location of the **ASSKEA proVisio^{Basic} or ASSKEA proVisio^{Basic} battery**. If the essential functions of the device are affected by electromagnetic interference, a reduced or non-existent volume flow and a reduced or non-existent vacuum build-up can be expected.

17.4 Overview of all cables and transducers that can be replaced by the user

Designation	Specification	Maximum length
Power supply unit incl. connection cable	Type: <u>GTM46402-3713.8</u> Technical data: In: AC 100 – 240 V ~ 50-60 Hz / 1,0 A Out: DC 13,8 V / 2,68 A / 37 W	1.20 m

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