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Copyright © 2019 ASSKEA GmbH, Gebesee.

The safety of the **ASSKEA secretion aspirators M- and S-series** complies with the acknowledged rules of technology and the requirements of the **German Medical Devices Act**.

The **ASSKEA secretion aspirators M- and S-series** bear the **CE marking CE1434** in accordance with EU Council Directive 93/42/EEC concerning medical devices and meet the essential requirements of Annex I of this directive.

The **ASSKEA secretion aspirators M- and S-series** have been tested in accordance with IEC 62353.

The **quality management system** applied by ASSKEA GmbH is certified in compliance with the relevant international standards.

The **ASSKEA secretion aspirators M- and S-series** are medical aspirators classified as class IIa in accordance with EU Council Directive 93/42/EEC, Annex IX.

Errors and omissions excepted.

Contents

1	User information	2	2.4	Warranty	15
1.1	Using this instruction for use	2	3	Operation	16
1.2	Icons	3	3.1	Set-up and start-up	16
1.2.1	Device, packaging and accessories	3	3.1.1	Connecting the ASSKEA secretion aspirators S-series	17
1.2.2	Operating elements	4	3.1.2	Connecting the ASSKEA secretion aspirators M-series	17
1.3	Symbol convention	4	3.2	Start-up	18
1.4	Glossary	4	3.2.1	Connecting the reusable secretion canister system	19
1.5	Purpose	6	3.2.2	Connecting the disposable secretion canister system	20
1.5.1	Intended use ASSKEA secretion aspirators S-series	6	3.2.3	Suction tube	22
1.5.2	Intended use ASSKEA secretion aspirators M-series	6	3.3	Operation of the ASSKEA secretion aspirators S-series	22
1.5.3	Essential functions	6	3.3.1	Setting the vacuum	22
1.5.4	Applied parts	6	3.3.2	Aspiration	22
1.5.5	Indications	6	3.4	Operation of the ASSKEA secretion aspirators M-series	23
1.5.6	Contraindications	6	3.4.1	Display and LED status indicators	23
1.5.7	Restrictions on use	6	3.4.2	Setting the vacuum	25
1.6	Basic safety instructions – ATTENTION!	6	3.4.3	Aspiration	26
1.7	User requirements	9	3.5	Replacement of the disposable liner ‘OneWay’	27
1.8	Information on product liability	9	3.6.	Decommissioning	28
1.9	Material compatibility	9	4	Maintenance	29
2	Product description	10	4.1	Cleaning and care in outpatient and inpatient care	29
2.1	Overview of the ME system	10	4.1.1	General information	29
2.1.1	ASSKEA secretion aspirators S-series	10	4.1.2	Cleaning and disinfecting the surface of the device	29
2.1.2	ASSKEA secretion aspirators M-series	10	4.1.3	Cleaning and disinfecting the reusable secretion canister	30
2.2	Product contents	11	4.1.4	Cleaning and disinfecting the connecting tube	30
2.3	Product properties	11	4.1.5	Cleaning and disinfecting the external canister ‘Bag’	31
2.3.1	ASSKEA secretion aspirators S-series	11	4.1.6	Cleaning/disposal of the suction tube	31
2.3.2	ASSKEA secretion aspirators M-series	12	4.1.7	Cleaning and disinfecting the rinsing bottle	31
2.3.3	Information on the double filter system	13	4.2	Cleaning and care in homecare use	31
2.3.4	Information on the secretion canister systems	13	4.2.1	General information	31
2.3.5	Information on the battery of the ASSKEA secretion aspirators M-series	14	4.2.2	Cleaning the surface of the device	32
			4.2.3	Hygienic cleaning of the reusable secretion canister	32

4.2.4	Hygienic cleaning of the connecting tube	33
4.2.5	Hygienic cleaning of the external canister 'Bag'	33
4.2.6	Hygienic cleaning of the suction tube	33
4.2.7	Hygienic cleaning of the rinsing bottle	33
4.3	Reuse of the device	33
4.4	Maintenance and service	33
4.5	Inspection of the ASSKEA secretion aspirators M- and S-series	34
5	Troubleshooting	35
5.1	ASSKEA secretion aspirators S-series	35
5.2	ASSKEA secretion aspirators M-series	36
6	Transport, storage and disposal	38
6.1	Decontamination prior to shipment	38
6.2	Storage	38
6.3	Disposal	38
7	Technical data	39
7.1	ASSKEA secretion aspirators S-series	39
7.2	ASSKEA secretion aspirators M-series	40
8	EMC information	41
8.1	Electromagnetic environment, in which the ASSKEA secretion aspirators M- and S-series may be operated	42
8.2	Dealing with electromagnetic interactions	43
8.3	List of all cables and transformers replaceable by the user	43
9	Order information	44
10	Publishing information	45

1 User information

1.1 Using this instruction for use

Please read this entire instruction for use before operating the **ASSKEA secretion aspirators M- and S-series** for the first time. Alongside trained medical professionals and instructed relatives of patients, it is also intended that patients use the devices themselves, providing that they have been trained accordingly. In this case, all functions can be used safely.

If support is required for startup, operation or maintenance, please contact ASSKEA GmbH (see section 10). Please also report unexpected operating performance or incidents to ASSKEA GmbH.


Please read the safety instructions (see section 1.6) to avoid hazards.

This instruction for use is a component of the **ASSKEA secretion aspirators M- and S-series**. Therefore please keep this instruction for use in an easily accessible place. Include this instruction for use when passing the device on to third parties.

1.2 Icons

1.2.1 Device, packaging and accessories

Symbol	Meaning
	Attention: potential for bodily injury, health risks or possible property damage
	Follow the instruction for use
	Power supply unit
	Protection class II
	Protection class: Type BF (Body Floating)
	Temperature limitation
	Air pressure limitation
	Humidity limitation
	Protect from moisture
	Do not use if packaging is damaged!
	Follow the instruction for use

Symbol	Meaning
	Note with useful information and tips
	This device must not be disposed of in domestic waste.
	Order number
	Serial number
	Lot number
	Date of manufacture
	Manufacturer
	Operating time (power-on time and power-off time)
	Do not reuse.
	Fuse
	IP22 (see section 1.4)

1.2.2 Operating elements

Symbol	Meaning
	Adjustment direction for vacuum setting using the vacuum regulator
	LED status indicator for external power supply (only M-series)
	power switch setting 'ON' (only M20 and S-series)
	LED status indicator for double filter system (only M20D and M30D)
	LED status indicator for the rechargeable battery (only M-series)
	LED status indicator ON/OFF (only M-series)
	power switch setting 'OFF' (only M20 and S-series)

1.3 Symbol convention

Symbol	Meaning
•	Enumeration
1. 2.	Perform the actions in the order indicated.

1.4 Glossary

A

approx. Abbreviation for 'approximately'

Aspirate Aspirate is the generic term for secretions, bodily fluids and flushing liquids that are typically accumulated during aspiration of the upper airways. These can be aspirated easily using the devices described here.

C

Contamination Contamination means that bacteria and viruses from the aspirate have come into contact with the device.

D

DFS® Double filter system
An external bacterial filter and a bacterial filter integrated into the aspirator make up the double filter system. The double filter system effectively protects the interior of the device from contamination and overflow. It enables safe processing and rapid reuse of the product.

Degree of protection The degree of protection indicates the protection of applied parts against electric shock. Applied parts of type BF must be installed isolated from earth and are not suitable for direct application to the heart.

I

incl. Abbreviation for 'including'

IP22 International Protection/Protection Class
The Protection Class defines the degree of protection of the device against contact and ingress of liquids.

The **ASSKEA secretion aspirators M- and S-series** are protected against finger access and falling water drops at an inclination of up to 15°.

M

ME system Abbreviation for 'medical electrical system'

MRI Abbreviation for 'magnetic resonance imaging'
This technique generates sectional images of the human body with the aid of a very strong magnetic field for the purpose of analyzing the organs.

O

Overflow Overflow means that aspirate is sucked into the interior of the device.

P

Processing The processing procedure is required for each new patient. The term processing means that all parts that have (potentially) come into contact with the aspirate are cleaned, disinfected and replaced, if necessary. The processing procedure must only be performed by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH.

R

resp. Abbreviation for 'respectively'

1.5 Purpose

1.5.1 Intended use ASSKEA secretion aspirators S-series

The **ASSKEA secretion aspirators S-series** are mains-operated medical aspirators used for temporary aspiration of aspirate from the trachea. The **ASSKEA secretion aspirators S-series** are ideally suited for use in hospitals, but can also be used for care and rehabilitation in the private sphere (homecare) for the respective patient.

1.5.2 Intended use ASSKEA secretion aspirators M-series

The **ASSKEA secretion aspirators M-series** are network-independent, mobile medical aspirators and are used for temporary and primarily spontaneous aspiration of aspirate from the trachea. Typical areas of application are:

- in outpatient and inpatient care,
- in homecare, especially for aspiration in tracheostomy patients.

1.5.3 Essential functions

- Generation of vacuum (high vacuum)
- Generation of flow (S20, S20K, M20, M20D: low flow; S30, M30D: high flow)

1.5.4 Applied parts

The suction tube is an applied part of type BF.

1.5.5 Indications

- Tracheostomy patients
- Aspiration on impaired respiratory function
- Aspiration of blood, secretion and food residues from the oral cavity, the pharyngeal zone and the bronchial system
- Evacuation of vacuum mattresses and vacuum splints

1.5.6 Contraindications

The **ASSKEA secretion aspirators M- and S-series** are contraindicated for the following applications:

- Liposuction
- Gynecology applications
- Thoracic drainage
- Continuous drainage
- Applications in wound areas
- Dental applications

1.5.7 Restrictions on use

- In medical rooms where potential equalization is required (e.g. heart surgery)
- In hazardous areas/in the MRI environment
- Outdoors/during transport

1.6 Basic safety instructions – ATTENTION!

Warning of damage due to improper power supply

Improper operation causes excess voltage in the device, which may be transmitted to the operator.

- Ensure prior to start-up that the mains supply for connecting the **ASSKEA secretion aspirators S-series** is suitable for 230V AC with a mains frequency of 50 Hz and for connecting the **ASSKEA secretion aspirators M-series** that it is suitable at 100 V to 240 V AC with a mains frequency of 50 - 60Hz.
- Ensure prior to start-up in UL listed markets such as the USA or Canada that the mains supply is designed for a supply voltage of 120V AC.
- Use only the supplied power supply unit (type: GTM91099-6015-3.0-T2) with the **ASSKEA secretion aspirators M-series** resp. the supplied mains cable (type: H05VV-F 2x1,0 mm²) with the **ASSKEA secretion aspirators S-series**.

Risks associated with use in conditions that are not approved

- The devices are not intended to be used outdoors/during transport.
- The devices are not intended to be used in medical rooms where potential equalization is required (e.g. heart surgery).

- The devices are not intended to be used in hazardous areas/in the MRI environment.

Health risks associated with handling infectious agents or pathogens

Infectious agents and pathogens in the aspirate cause health risks.

- Never operate the device without an external bacterial filter.
(exception: ASSKEA disposable secretion canister system)
- Check that the external bacterial filter is dry and clean prior to each use to ensure faultless function! (exception: ASSKEA disposable secretion canister system)
- Always perform the aspiration using a suitable sterile disposable catheter. The suction tube must never come into direct contact with the aspiration site.
- Follow the hygiene, cleaning and decontamination instructions (see section 4).

Risk of strangulation

- People may strangle themselves on the tubing or the mains cable, especially if tubes resp. cables are unduly long.
- Ensure that no unauthorized/uninvolved persons are near the device during aspiration.
- Store the device incl. accessories in the shipping carton until the next time of use.

Damage to the device due to heat build up

- Do not cover the power supply unit.
- Keep the aspirator as well as the mains cable and the power supply unit away from other heat sources.
- Do not place the aspirator directly adjacent to other devices, since this can result in significant warming of the aspirator.

Safety defects due to incorrect accessories and spare parts

The use of accessories and spare parts, detachable parts or materials that are not recommended by ASSKEA GmbH and specified in the instruction for use may compromise the safety and function of the device. Any warranty is excluded for damage caused by using accessories and spare parts not recommended or by improper use. Use only the recommended original accessories and spare parts.

Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Prior to each use, check the housing for any damage and do not operate the device if there is obvious damage to the housing.
- Prior to each use, check the suction tube, secretion canister, bacterial filter and, if ap-

plicable, any other accessories that are subject to wear and damage, to ensure that the components are in perfect condition and proper operation of the device is guaranteed. If this is not the case, replace the parts immediately.

Risks to people due to improper handling

- Use the device only for its intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. thoracic drainage).
- If aspiration procedures are carried out too often, slight bleeding may occur.
- When using the power supply unit, ensure that it is connected to the aspirator first and then to the mains supply (100V to 240V AC). Disconnection of the power supply unit from the mains supply must be performed in exactly the opposite sequence (first remove the power supply unit from the mains supply (100 V to 240 V AC) and then remove it from the aspirator).
- When using the mains cable, ensure that it is connected to the aspirator first and then to the mains supply (230 V). Disconnection of the mains cable must be performed in the opposite sequence.
- Never touch parts of non-ME devices in the patient environment and the patient simultaneously.

Risk of damage due to electromagnetic phenomena

Medical electrical devices are subject to special precautionary measures regarding electromagnetic compatibility and must be installed and operated in accordance with the EMC information provided in the accompanying documents (see section 8)!

Warning of safety defects due to improper connections of the ME system

Connection of the ME system with other devices or installations or pieces of equipment not recommended by ASSKEA GmbH and not specified in the instruction for use may compromise the safety and function of the ME system. Any warranty is excluded for damage caused by connecting devices or installations or pieces of equipment not recommended with the **ASSKEA secretion aspirators M- and S-series** or by improper use. Connect only recommended original parts with the **ASSKEA secretion aspirators M- and S-series**.

Damage to the device from ingress of liquids

The **ASSKEA secretion aspirators M- and S-series** have the IP classification IP22 against the ingress of liquids. Nevertheless protect the device from moisture (see symbol in section 1.2.1).

- Do not use the device near splashing water.
- Do not use the device in damp rooms or

while bathing or showering.

- Do not allow the mains cable or power supply unit, power switch, socket for the power supply unit/mains cable or the display foil to get wet.
- Never submerge the device in water or other liquids (even when it is switched off).

Inspection of the internal power supply

Since the internal battery of the **ASSKEA secretion aspirators M-series** is not automatically kept in a fully operational condition, the state of charge must be checked periodically and, if necessary, the battery must be replaced by service staff. The battery may be replaced only by authorized service staff, since replacement by inadequately trained persons could result in a hazard (such as excessive temperatures, fire or explosion)!

Known or identifiable conditions when using the device for medical care within a home environment

- Children and pets must be kept away from the device. Ensure that the device cannot be knocked over or fall in such a way that persons can be hit by the falling device.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the home's power supply.
- Do not use the device in damp rooms, bathrooms or showers.

- Do not allow the mains cable or power supply unit, power switch, socket for the power supply unit/mains cable or the display foil to get wet.
- Never submerge the device in water or other liquids (even when it is switched off).
- Ensure that lint and dust are immediately removed from the device and accessories in order to ensure unimpaired functionality. In addition, vermin must not be allowed in the proximity of the device, as there is a risk they may get inside the device and cause damage.
- Never expose the device or accessories to direct sunlight as this could result in a strong heating effect, causing functional impairments.
- Some devices and sources that are commonly used in the home environment may cause potential faults in the device or accessories, e.g. fireplaces or radiant heaters (strong heating of the device) or inhalers or steam kettles (excessive air humidity). Do not operate such devices and sources in the vicinity of the **ASSKEA secretion aspirators M- and S-series**.

Possible physiological effects and unobvious risks

- Select the vacuum depending on the respective patient and on the medical indication to avoid personal injuries.

- Never use the **ASSKEA secretion aspirators M- and S-series** to aspirate solids or inhaled foreign objects out of the throat. This could lead to the object being pushed deeper into the airways potentially causing total obstruction.
- Always place the device upright on a sturdy and flat base without slope. Make sure that the device cannot be knocked over or fall in such a way that persons could be hit by the falling device.
- Other devices, examinations or treatments may be influenced by the device. For this reason, special attention should always be paid to other devices as well as to examinations and treatments performed in parallel so that any effect is detected as quickly as possible.
- When using the device, adequate lighting must always be ensured so that all labels can be identified clearly.
- Small, detached parts could be inhaled or swallowed. Therefore ensure that no unauthorized persons, children or pets are near the device.
- Although the materials used have been tested for compatibility, in exceptional cases allergic reactions to the exposed materials on the device and the accessories may occur. This applies especially for contact injuries after prolonged exposure. In such cases, consult a doctor immediately.

1.7 User requirements

The **ASSKEA secretion aspirators M- and S-series** must only be operated and used by instructed and appropriately trained persons. Familiarize yourself with the mode of operation of the **ASSKEA secretion aspirators M- and S-series** prior to start-up.

Training on how to use the **ASSKEA secretion aspirators M- and S-series** is provided by ASSKEA GmbH or an authorized distribution partner of ASSKEA GmbH. Product training takes approx. one to two hours and, among other things, includes an explanation of the design and function of the device, handling of the device, cleaning and disinfection as well as the procedure to be followed for each new patient and for disposal.

Training should be repeated regularly every 24 months.

Each participant receives a certificate as proof of training.

1.8 Information on product liability

Liability for the operation of the devices is transferred to the operator in the following cases:

- if the device is used for a purpose other than its intended use,
- if the device is not used in accordance with the instruction for use,
- if the device is opened by unauthorized personnel,
- if installation, settings, enhancements, maintenance or repairs are performed by unauthorized personnel,
- if incorrect accessories or spare parts are used,
- if the safety seal is removed/broken.

Advice for the responsible organization:

The assembly of ME systems and modifications during their actual service life require that observance of requirements of the applicable standards is regularly reviewed.

1.9 Material compatibility


Aggressive substances may damage the device and the accessories.

Please follow the cleaning and care instructions (see section 4).

2 Product description

2.1 Overview of the ME system


2.1.1 ASSKEA secretion aspirators S-series (S20K, S20, S30)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see figs. 29 and 30)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Vacuum meter
- E** Vacuum regulator
- F** Power switch
- G** Mains connection 230 V AC
- H** Mains cable
- I** Suction tube

Fig. 1


2.1.2 ASSKEA secretion aspirators M-series a) analog (M20)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see figs. 29 and 30)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Vacuum meter
- E** Vacuum regulator
- F** Power switch
- G** Socket for power supply unit (12 V DC)
- H** Power supply unit (GTM91099-6015-3.0-T2)
- I** Mains cable
- J** Suction tube

Fig. 2

b) digital (M20D, M30D)

- A** Secretion canister (reusable) with integrated overflow protection in the lid (see figs. 29 and 30)
- B** Bacterial filter (external)
- C** Connecting tube
- D** Buttons for vacuum adjustment
- E** ON/OFF button
- F** Max. button
- G** Status indicator
- H** Socket for power supply unit (12 V DC)
- I** Power supply unit (GTM91099-6015-3.0-T2)
- J** Mains cable
- K** Suction tube

Fig. 3

2.2 Product contents

- **ASSKEA secretion aspirator S-series** or **ASSKEA secretion aspirator M-series**
- instruction for use
- mains cable (for **ASSKEA secretion aspirators S-series**)
- the power supply unit GTM91099-6015-3.0-T2 incl. mains cable (for **ASSKEA secretion aspirators M-series**)
- reusable secretion canister system or disposable secretion canister system, if applicable
- charging instructions, multilingual (only for **ASSKEA secretion aspirators M-series**)
- instructions for safe handling of battery packs (only for **ASSKEA secretion aspirators M-series**)
- 'Used Medical Device' label and decontamination certificate
- test report according to IEC 62353
- optional accessories (depending on the order)

2.3 Product properties

Hazard to persons due to improper handling

- Use the device only for the intended purpose.
- Never use the device for aspiration in the low vacuum range (e.g. thoracic drainage).

Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Do not operate the device if there is obvious damage to the housing.

2.3.1 ASSKEA secretion aspirators S-series

The **ASSKEA secretion aspirators S-series** are exceptionally convenient aspiration devices for stationary use.

The devices are operated via a mains cable to the 230 V mains power supply. The desired vacuum is generated by a maintenance-free vacuum pump. An overtemperature protection prevents the aspiration device from overheating and leads to its automatic switch-off. After it is switched on via the power switch, the vacuum pump generates a vacuum

in the tubing system and in the secretion canister, which is used to aspirate secretion (via a suitable sterile aspiration catheter). The aspirate is directed away from the patient and collected in the secretion canister. Once the secretion canister is full, aspiration is stopped by means of a mechanical overflow protection integrated into the lid of the reusable secretion canister or a self-sealing bacterial filter in the disposable liner 'OneWay'.

Flow rate is adjusted by using the vacuum regulator.

The reusable secretion canister or holder for the external canister 'Bag' is attached by directly pushing it onto the guide rails on the housing of the aspiration device (figs. 11, 12). The vacuum connection between the aspiration device and the secretion canister is made via the connecting tube (figs. 1, 2, 3(C)).

2.3.2 ASSKEA secretion aspirators M-series

The **ASSKEA secretion aspirators M-series** are exceptionally convenient, battery-powered aspiration devices for mobile use.

All mobile devices are operated via the internal battery or via the supplied power supply unit that also can be used to recharge the battery. The desired vacuum is generated by a maintenance-free vacuum pump. An over-temperature protection prevents the battery from overheating during operation and while charging under unfavorable environmental conditions.

After it is switched on, the vacuum pump generates a vacuum in the tubing system and in the secretion canister, which is used to aspirate secretion (via a suitable sterile aspiration catheter). The aspirate is directed away from the patient and collected in the secretion canister. Once the secretion canister is full, aspiration is stopped by means of a mechanical overflow protection integrated into the lid of the reusable secretion canister or a self-sealing bacterial filter in the disposable liner 'OneWay'.

The reusable secretion canister or holder for the external canister 'Bag' is attached by directly pushing it onto the guide rails on the housing of the aspiration device (figs. 11, 12). The vacuum connection between the aspira-

tion device and the secretion canister is made via the connecting tube (figs. 1, 2, 3(C)).

The **ASSKEA secretion aspirators M-series** are fundamentally divided according to the configuration versions into the following groups:

a) analog (M20 - fig. 4)

Similar to stationary aspiration devices, the vacuum is adjusted by using the vacuum regulator (additional air regulation) in the analog, mobile version. All **ASSKEA secretion aspirators M-series** are operated via the internal battery or via the supplied power supply unit that also can be used to recharge the battery.

b) digital (M20D, M30D - fig. 5)

The digital, mobile versions are equipped with a keyboard foil that can be used to adjust the desired vacuum value and to regulate it using an electronic sensor. All **ASSKEA secretion aspirators M-series** are operated via the internal battery or via the supplied power supply unit that also can be used to recharge the battery.


Fig. 4


Fig. 5

2.3.3 Information on the double filter system

Health risks associated with handling of infectious agents or pathogens

Infectious agents and pathogens in the aspirate cause health risks.

- Never operate the device without an external bacterial filter.
(exception: ASSKEA disposable secretion canister system)
- Check that the external bacterial filter is dry and clean prior to each use to ensure faultless function!
(exception: ASSKEA disposable secretion canister system)

The ASSKEA double filter system DFS® consists of the external bacterial filter and the bacterial filter integrated into the aspirator.

For the reusable secretion canister system, a separate, external hydrophobic (liquid-repellent) bacterial filter is used.

For the disposable secretion canister system, an external self-sealing bacterial filter is integrated into the disposable liner 'OneWay'. The external bacterial filters (reusable as well as disposable secretion canister system) and the integrated bacterial filter form the ASSKEA double filter system and in combination they are 99.999% effective against bacteria and viruses.

The ASSKEA double filter system DFS® effectively protects against overflow and contamination of the interior of the device. It permits fast, simple and cost-effective processing.

Service life and reuse

The **external bacterial filter** (for the reusable secretion canister system) is intended for **single use**. For this reason, it must be replaced for each new patient. When used exclusively for one patient, replace the external bacterial filter no less than every two weeks, if it is used frequently.

When replacing the external bacterial filter, follow the supplied instruction for use!

The **integrated bacterial filter is not suitable for reuse**. To ensure consistent performance, the integrated bacterial filter must be replaced by authorized service staff **prior to each new patient**, if it has come into **contact with the aspirate** (blocking) or during **maintenance/repair**. This also applies after the filter service life (see below) has expired.

For the digital versions **M20D** and **M30D** of the **ASSKEA secretion aspirators** an **LED** additionally helps to indicate the filter status of the double filter system (DFS).

If the **DFS LED** illuminates green, the standby function of the filter system is guaranteed. If it illuminates red, the integrated or external bacterial filter may be blocked or aspiration is hindered by blockage in the tubing system. Before replacing the **DFS**, first check the aspiration system for possible blockages.

If the integrated bacterial filter is blocked, it must be replaced by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH.

2.3.4 Information on the secretion canister systems

Disposable secretion canister system:

The disposable secretion canister system consists of the external canister 'Bag', the holder for the external canister 'Bag', the disposable liner 'OneWay', the connecting tube and the suction tube. The disposable liner 'OneWay' has an integrated bacterial filter, carbon filter and solidifier. The self-sealing bacterial filter integrated in the disposable liner 'OneWay' is 99 % effective against microorganisms. It prevents overflow in the event of an operational error. If liquid reach-

es this filter, aspiration is no longer possible and the aspiration process is discontinued. The disposable liner 'OneWay' must be replaced.

The activated carbon filter in the disposable liner 'OneWay' reduces the spread of odor.

Solidifier

Disposable liners 'OneWay' filled with aspirate can be transported and disposed of in a leak-proof manner by using the solidifier. The aspirate solidifies after an average gelling time of 2-5 minutes (depending on the consistency of the secretion), regardless of the aspiration intervals.

The **disposable liner 'OneWay' and the suction tube** are intended for **single use**. Replace the disposable liner 'OneWay' and suction tube in accordance with the respectively applicable hygiene instructions, but no less than once per week, prior to each new patient or if it is full.

Reusable secretion canister system:

The reusable secretion canister system consists of the reusable secretion canister, the lid of the secretion canister with integrated overflow protection, the suction tube and a connecting tube with external bacterial filter.

The external bacterial filter is 99.999% effective against bacteria and viruses. It prevents an overflow in the event of operational error. If liquid reaches this filter, aspiration is no longer possible and the aspiration process is discontinued. The external bacterial filter must be replaced.

2.3.5 Information on the battery of the ASSKEA secretion aspirators M-series

Inspection of the internal power supply

- Since the internal battery is not automatically kept at a fully operational condition, the state of charge must be checked periodically and, if necessary, the battery must be replaced by service staff.
- The battery may be replaced only by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH, since replacement by inadequately trained persons could result in a hazard (such as excessive temperatures, fire or explosion)!

The state of charge of the battery is indicated by the corresponding LED status indicator (see section 1.2.2).

It is strongly recommended to fully charge the battery prior to the first start-up of the **ASSKEA secretion aspirators M-series** and to repeat this after the first uses. The typical service life of the battery is approx. 500 load cycles; approx. 70 % of the primary capacity is still available thereafter.

The **ASSKEA secretion aspirators M-series** are equipped with a lithium ion battery, which unlike traditional types of rechargeable batteries has a low self-discharge rate. In addition lithium ion batteries do not have a memory effect. Therefore they can and should be recharged at any time after initial charging. Only frequent short-time charging should be avoided.

The **ASSKEA secretion aspirators M-series** should ideally be stored and charged at room temperature in accordance with the ambient conditions specified in the technical data. Never store the devices and battery in a discharged state!

Fully recharge the battery, if the devices are not operated for a prolonged period of time (approx. 10 months).

The battery of the **ASSKEA secretion aspirators M-series** is protected against depth discharge by means of protective measures; nevertheless the charging information listed

above must be followed. The battery is also protected against overheating during charging. If the battery temperature is exceeded during charging due to improper ambient conditions, charging is temporarily discontinued to allow cooling. The purpose of this measure is to ensure safe operation and to protect the battery.

2.4 Warranty

WARNING: Modification of the ME device is not permitted!

The **ASSKEA secretion aspirators M- and S-series** are covered by warranty for 2 years. It is neither extended nor renewed by warranty work. The battery is covered by warranty for 6 months. Wearing parts are excluded from the warranty.

ASSKEA GmbH is responsible for impacts on safety, reliability and specified performance, only if:

- original ASSKEA accessories and spare parts are used,
- maintenance and repair are performed by professionals authorized by ASSKEA GmbH or by ASSKEA itself,
- the affected product is used and operated

in accordance with the instruction for use and within its intended use.

All warranty claims are voided, if:

- the device is opened by unauthorized persons,
- the safety seal is removed/broken,
- repairs are performed by unauthorized persons,
- modifications are made to the device,

since the basic safety of the device can no longer be guaranteed and functional limitations may occur in these cases.

3 Operation

Hazard to persons due to improper handling

- Please read through sections 3.1 and 3.2!
- The **ASSKEA secretion aspirators M- and S-series** may be operated and used only by instructed and appropriately trained persons.
- Use only suitable sterile disposable catheters for aspiration! Commercially available sterile disposable catheters from CH10 to CH16 that are intended for the use in the tracheal area can be connected.
- If aspiration procedures are carried out too often, slight bleeding may occur.

Malfunction due to aspirated secretion

- Make sure that you empty the reusable secretion canister when it is half-full to curtail frothing. If the reusable secretion canister is full, the mechanical overflow protection is triggered. This disrupts the aspiration process.
- Ensure that the disposable liner 'OneWay' is replaced in a timely manner. The integrated overflow

protection is triggered if the disposable liner 'OneWay' is full. This disrupts the aspiration process.

- Always have one backup disposable liner 'OneWay' ready or when using the reusable secretion canister system one backup external bacterial filter, since these are absolutely necessary for safe operation!
- Switch off the device when emptying the reusable secretion canister or when replacing the disposable liner 'OneWay'.
- If secretion has been aspirated into the device, it must be properly processed by ASSKEA GmbH or by a service partner authorized by ASSKEA GmbH!

Damage to the devices due to insufficient acclimatization

After the devices have been exposed to temperatures between -20 °C and +60 °C during transport/storage according to the technical data (see section 7), they must first acclimatize for ca. 2 h at room temperature (ca. 20 °C) before the use as intended is possible.

3.1 Set-up and start-up

Known or identifiable conditions when using the device for medical care within a home environment

- Children and pets must be kept away from the device to ensure that it cannot be knocked over or dropped.
- Prior to connecting the power supply unit, ensure that the voltage of the device corresponds to the home's power supply.
- Do not use the device in damp rooms, bathrooms or showers.
- Do not allow the mains cable or power supply unit, power switch, socket for the power supply unit/ mains cable as well as the display foil to get wet.
- Never submerge the device in water or other liquids (even when it is switched off).

The following sections describe the operating elements, connections and the start-up of the **ASSKEA secretion aspirators M- and S-series**:

3.1.1 Connecting the ASSKEA secretion aspirators S-series

Check whether the stated mains cable (see section 7 'Technical data') is present and whether there is any damage. Replace the mains cable immediately if there is any damage.

For the mains-operated aspiration devices, use the mains connection (fig. 1(G)) to connect the device to the mains power supply with the aid of the mains cable. Make sure that the device is positioned in such a way that it is possible to easily disconnect it later. First connect the mains cable to the mains connection of the **ASSKEA secretion aspirator S-series** and then to the mains power supply. Information on the permissible environmental conditions during operation can be found in section 7 'Technical data'.


Fig. 6 Connection of the ASSKEA secretion aspirators S-series to patient and accessories

3.1.2 Connecting the ASSKEA secretion aspirators M-series

Check whether the stated power supply unit and mains cable (see section 7 'Technical data') is present and whether there is any damage. Replace the power supply unit or mains cable immediately if there is any damage.

For the network-independent aspiration devices, use the socket for the 12 V power supply unit (fig. 2(G), fig. 3(H)) to connect the device to the mains power supply with the aid of the supplied power supply unit to charge or operate it as needed. Make sure

that the device is positioned in such a way that it is possible to easily disconnect it later. First connect the power supply unit to the socket for the power supply unit of the **ASSKEA secretion aspirator M-series** and then to the mains power supply. Information on the permissible environmental conditions during operation can be found in section 7 'Technical data'.


Fig. 7 Connection of the ASSKEA secretion aspirators M-series to patient and accessories

3.2 Start-up

The following explanations for start-up apply to all versions of the **ASSKEA secretion aspirators M- and S-series**.

- Please read this entire instruction for use prior to initial start-up of the **ASSKEA secretion aspirator M- and S-series**.
- Remove the device and the accessories from the packaging.
- Always place the device on a sturdy and flat surface without slope.
- The **ASSKEA secretion aspirators M- and S-series** can be placed next to the patient bed. Optionally, a variable holder for the attachment of the devices to tube and rail systems as well as a bed holder are available.
- Inspect all tubes as well as the power supply unit or mains cable for damage prior to each start-up of the **ASSKEA secretion aspirators M- and S-series**. It is important to avoid kinking when connecting the tubes. Ensure prior to switching on the unit that the secretion canister and the tubing are properly connected.
- Fully charge the battery prior to initial start-up (only valid for **ASSKEA secretion aspirators M-series**).


- Prior to the initial start-up it is essential that you follow the safety instructions in section 1.6.
- Always have one backup disposable liner 'OneWay' ready when using the reusable secretion canister system or one backup external bacterial filter, since these are absolutely necessary for safe operation!

To ensure optimum aspiration, place the **ASSKEA secretion aspirators M- and S-series** below the patient to be treated. It is important that the suction tube does not form a dip and is situated at least at patient level.


3.2.1 Connecting the reusable secretion canister system

1. Screw the lid of the reusable secretion canister sufficiently tight and attach the canister to the device.


2. Insert the connecting tube with the blue elbow fitting toward the side of the device. Make sure that the side of the external bacterial filter with the lettering 'IN' is facing the canister (directed toward patient)!

Fig. 8

3. Place the grey elbow adapter of the connecting tube onto the connection in the middle of the reusable secretion canister lid (above the overflow protection).

Fig. 9

4. Place the suction tube onto the second connection (patient connection). Ensure that the tube is fixed sufficiently tight.

Fig. 10

The reusable secretion canister or the holder for the external canister 'Bag' is attached by directly pushing it onto the guide rails on the housing of the device (figs. 11, 12).

Please note that up to four mounting positions can be used for the secretion canister (figs. 11, 12) on all devices.


Fig. 11

Fig. 12

3.2.2 Connecting the disposable secretion canister system



Malfunction due to collapsing disposable liner 'OneWay'

A leak in the external canister 'Bag' or at the lid of the disposable liner 'OneWay' resp. damage of the canister or the T-piece may cause air to flow into the external canister 'Bag'. This may lead to the collapse of the disposable liner 'OneWay'.

- Inspect the disposable secretion canister system to ensure that the lid of the disposable liner 'OneWay' is securely connected to the external canister 'Bag'.
- Ensure that all connections are firmly attached and properly connected.
- Ensure that the external canister 'Bag' is undamaged and the T-piece is securely attached.
- Follow the instruction for use supplied by the manufacturer!

The original ASSKEA disposable secretion canister system consists of the external canister 'Bag', the holder for the external canister 'Bag', the disposable liner 'OneWay', the connecting tube and the suction tube.


Fig. 15

- Place the disposable liner 'OneWay' in the reusable external canister 'Bag'. Press the lid's edges down firmly to ensure proper sealing.


Fig. 17

- Connect the vacuum connection of the device with the corresponding vacuum connection of the external canister 'Bag' (top end of the T-piece). Use the supplied connecting tube to do so.


Fig. 16

- Attach the prefitted connecting tube of the disposable liner 'OneWay' to the bottom end of the T-piece located at the external canister 'Bag'.


Fig. 18

- Connect the patient connection of the disposable liner 'OneWay' (fig. 13(B)) to the suction tube. Ensure that the tube is attached sufficiently securely.

3.2.3 Suction tube

For all of the **ASSKEA secretion aspirators M- and S-series** you can connect the suction tube with the aspiration catheter using a fingertip.

3.3 Operation of the ASSKEA secretion aspirators S-series

3.3.1 Setting the vacuum


Fig. 19

1. Set the power switch (fig. 1(F)) to 'I'.
2. Keep the opening of the suction tube (fig. 19) and the fingertip (fig. 20) shut.
3. Turn the vacuum regulator (fig. 1(E)) until the vacuum meter (fig. 1(D)) displays the desired vacuum. Please note that higher vacuum values are set by


Fig. 20

- turning the vacuum regulator clockwise and that the indicator of the manometer moves counter clockwise at the same time.
4. You can also perform fine adjustments to the set vacuum via the additional air opening on the fingertip.

3.3.2 Aspiration

1. Perform the aspiration procedure as you were instructed by specialized personnel.
2. Control the aspiration procedure and the desired vacuum using the additional air opening on the fingertip (fig. 20).
3. Shut down the device after the aspiration procedure according to section 3.6.

Please note that the pump can not be started beyond an existing vacuum of -0.5 bar when switching the device on again.

Please note that the **ASSKEA secretion aspirators S-series** are designed for **short-term** operation (20 minutes)! Please also note the technical data for this.

Always take note of additional important instructions in section 2.3.3 regarding handling of the ASSKEA double filter system!

3.4 Operation of the ASSKEA secretion aspirators M-series

3.4.1 Display and LED status indicators

a) analog secretion aspirator M20


Fig. 21

Symbol	Description	Meaning of display
	LED illuminates continuously green	Power supply unit connected; voltage available
	LED illuminates continuously yellow	Power supply unit connected; power supply permissible, but not corresponding to the intended input voltage of 12V DC → Battery is being charged!
	LED illuminates continuously red	Input voltage outside the permissible range; automatic switch-off to protect the device
	LED illuminates continuously green	Battery capacity 100% to 40%
	LED illuminates continuously yellow	Battery capacity < 40% to 10%
	LED illuminates continuously red	Battery capacity < 10% to 0%, immediately plug in the external power supply. After completing the aspiration procedure, turn off the device for charging!
	LED illuminates continuously red in spite of external power supply	Battery faulty; operation possible only using the plugged in power supply unit
	LED blinks green	Battery is being charged.
	LED blinks red	Overtemperature protection; interruption of the charging/discharging process until the temperature of the battery is within permissible range again
	LED illuminates continuously green	Pump is switched on.

b) digital secretion aspirators M20D and M30D


Fig. 22

Symbol	Description	Meaning of display
	LED illuminates continuously green	Power supply unit connected; voltage available
	LED illuminates continuously yellow	Power supply unit is connected; power supply permissible, but not corresponding to the intended input voltage of 12V DC → Battery is being charged!
	LED illuminates continuously red	Input voltage outside the permissible range; automatic switch-off to protect the device
	LED illuminates continuously green	Battery capacity 100% to 40%
	LED illuminates continuously yellow	Battery capacity < 40% to 10%
	LED illuminates continuously red	Battery capacity < 10% to 0%, immediately plug in the external power supply. After completing the aspiration procedure, turn off the device for charging!
	LED illuminates continuously red in spite of external power supply	Battery faulty; operation possible only using the plugged in power supply unit
	LED blinks green	Battery is being charged.
	LED blinks red	Overtemperature protection; interruption of the charging/discharging process until the temperature of the battery is within permissible range again
	LED illuminates continuously green	Double filter system operational
DFS	LED illuminates continuously yellow	Filter monitoring time has been reached. The integrated bacterial filter should be replaced and the monitoring time should be reset promptly. For this purpose, return the device to ASSKEA GmbH or to one of the service partners authorized by ASSKEA GmbH.
	LED illuminates continuously red	DFS [®] blocked; the double filter system must be replaced; in particular, the integrated bacterial filter must be replaced. Adherence detection: the tubing system must be checked for blockages.

Symbol	Description	Meaning of the display
	LED illuminates continuously green	Pump is switched on.
LED -0.1 bar to	LED illuminates red for 2 seconds	Display of the set or preset vacuum (target value display); factory-provided default: -0.4 bar
-0.8 bar	LED illuminates continuously green	Display of the current vacuum (actual value display)

3.4.2 Setting the vacuum

a) analog secretion aspirator M20


Fig. 23


Fig. 24

1. Set the power switch (fig. 2(F)) to 'I'.
2. Keep the opening of the suction tube (fig. 23) and of the fingertip (fig. 24) shut.
3. Turn the vacuum regulator (fig. 2(E)) until the vacuum meter (fig. 2(D)) displays the desired vacuum. Please note that higher vacuum values are set by turning the vacuum regula-

- tor clockwise and that the indicator of the manometer moves counter clockwise at the same time.
4. You can also perform fine adjustments to the set vacuum using the additional air opening on the fingertip.

b) digital secretion aspirators M20D and M30D

The digital network-independent versions operate vacuum-controlled. The vacuum is set via the keyboard foil on the front side of the devices, where the tolerance for the accuracy of the vacuum display is 10%. This technology enables efficient and economical operation with the energy resources during battery operation.


Fig. 25

1. Switch on the device by pressing the ON/OFF button (fig. 3(E)). The ON LED signalizes that the pump is switched on by illuminating green.

Immediately after switching on, one of the 8 vacuum display LEDs (fig. 3 (G)) illuminates red for 2 seconds to show the vacuum value that was last set or that was present during start-up (target value display).

2. Set your desired vacuum via the arrow keys (fig. 26(1)). When the device is later switched off, the last set value is saved.


Fig. 26

3. The display via the red LED signalizes the selected target value of the vacuum to you. The automatic switch over to the current vacuum via a green LED occurs after 2 seconds.

4. You can also perform fine adjustments to the set vacuum via the additional air opening on the fingertip.


Fig. 27

For immediate maximum vacuum:

- If you require the maximum value quickly, press the 'max' button (fig. 28(2)). You can also set the maximum vacuum by pressing the 'up' button several times.
- Preselection of the maximum vacuum is signalized via a red LED for 2 seconds (fig. 28(3)).


Fig. 28

3.4.3 Aspiration

1. Perform the aspiration procedure as you were instructed by specialized personnel.
2. Control the aspiration procedure and the desired vacuum using the additional air opening on the fingertip (fig. 27).
3. Shut down the device after the aspiration procedure according to section 3.6.

Please note that, when using the analog M20, the pump can not be started beyond an existing vacuum of -0.5 bar when switching the device on again.

Always take note of additional important instructions in section 2.3.3 regarding handling of the ASSKEA double filter system!

Ensure that the battery has a sufficient capacity during the aspiration procedure. The battery status display must signal green or yellow. If the display turns red, the external power supply unit must be plugged in immediately. In this case, aspiration can continue without interruption.

3.5 Replacement of the disposable liner 'OneWay'

Health risks associated with handling infectious agents or pathogens

Infectious agents and pathogens in the aspirate cause health risks.

- Wear suitable disposable gloves when replacing the disposable liner 'OneWay'.
- Use the disposable liner 'OneWay' for one patient only.

- Replace the disposable liner 'OneWay' and suction tube in accordance with the respectively applicable hygiene instructions, no less than once per week, prior to each new patient or if it is full.
- For each new patient, professional hygienic processing by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH is strictly required!
- After each aspiration procedure, all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- Disposal of aspirate and contaminated components must be performed in a proper manner.

1. Switch off the device of the **ASSKEA secretion aspirators M- and S-series** by setting the power switch (fig. 1(F), fig. 2(F)) of the analog secretion aspirators to 'O' or by pressing the ON/OFF button (fig. 3(E)) of the digital secretion aspirators for 3 seconds (confirmed by an acoustic signal).
2. Separate the suction tube from the patient connection.
3. Separate the pre-assembled connecting tube of the disposable liner 'OneWay' at the bottom end of the T-piece of the external canister 'Bag'.

4. Remove the disposable liner 'OneWay' from the reusable external canister 'Bag'.
5. Dispose of the disposable liner 'OneWay' and suction tube in a proper manner (see section 6.3).
6. Insert a new disposable liner 'OneWay' in the reusable external canister 'Bag' as specified in section 3.2.2. Ensure that the connecting tube as well as the lid of the disposable liner are properly seated on the external canister.
7. Attach a new suction tube to the patient connection of the disposable liner 'OneWay'.
8. Switch on the **ASSKEA secretion aspirator M- and S-series**.

If the reusable secretion canister system is used, please follow the instructions regarding cleaning/disinfection in section 4.1.3 or 4.2.3 and 4.1.4 or 4.2.4.

3.6 Decommissioning

Health risks associated with handling infectious agents or pathogens

Infectious agents and pathogens in the aspirate cause health risks.

- Wear appropriate disposable gloves during decommissioning.
- Reuse of bacterial filters for multiple patients is prohibited for safety reasons!
- For each new patient professional hygienic processing by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH is strictly required!
- After each aspiration procedure all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- Disposal of aspirate and contaminated components must be performed in a proper manner.

1. Switch off the device after the aspiration by setting the power switch (fig. 1(F), fig. 2(F)) of the analog secretion aspirators to 'O' or by pressing the ON/OFF button (fig. 3(E)) of the digital secretion aspirators for 3 seconds (confirmed by an acoustic signal).
2. Disconnect the power supply unit or the mains cable from the mains supply (100V to 240V AC) and then disconnect the device plug from the aspiration device.
3. Remove the suction tube (fig. 1(I), 2(J), 3(K)) from the aspiration catheter and from the secretion canister and clean (see section 4.1.6 or 4.2.6) or dispose it of (see section 6.3) properly.
4. Remove the disposable liner 'OneWay' as described in section 3.5 up to and including point 5. Then remove the external canister 'Bag' from the belonging holder. Pull the holder for the external canister 'Bag' off the guide rails on the device and follow the instructions for cleaning the external canister 'Bag' in section 4.1.5 or 4.2.5.
If you are using the reusable secretion canister system, pull the reusable secretion canister off the guide rails on the de-

vice and then follow the instructions for cleaning or disposing of the components of the reusable secretion canister system described in section 4.1.3 or 4.2.3 and 4.1.4 or 4.2.4.

5. Clean the surface of the device according to section 4.1.2 or 4.2.2.
6. Store the device in the shipping carton until next use.

4 Maintenance

4.1 Cleaning and care in outpatient and inpatient care

4.1.1 General information

Health risks associated with handling infectious agents or pathogens

Infectious agents and pathogens in the aspirate cause health risks.

- Never operate the device without an external bacterial filter.
(exception: ASSKEA disposable secretion canister system)
- Replace the external bacterial filter prior to each new patient or at least every two weeks when used exclusively for one patient, if it is used frequently!
(exception: ASSKEA disposable secretion canister system)
- Reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
- Wear appropriate disposable gloves when replacing components that have come into contact with the aspirate.
- Replace the disposable liner 'One-Way' and suction tube in accordance with the respectively applicable hy-

giene instructions, no less than once per week, prior to each new patient or if it is full.

- For each new patient, professional hygienic processing by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH is strictly required! Follow the information in section 4.3.
- After each aspiration procedure all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- Disposal of aspirate and contaminated components must be performed in a proper manner.

Health risks associated with handling disinfectants

- Use of appropriate protective clothing during disinfection is recommended.
- Follow the instructions of the manufacturer of the disinfectant.

Possible bodily injury due to electric shock

- Switch off the device prior to cleaning/disinfection.
- Disconnect the power supply unit resp. the mains cable by unplugging it from the mains supply. Then dis-

connect the power supply unit or the mains cable from the socket for power supply of the **ASSKEA secretion aspirator M- and S-series**.

Damage to the device due to improper cleaning agents

- Do not use disinfectants that contain acetone. These may damage or disfigure the housing components and the accessories.
- Follow the instructions for use provided by the manufacturers of the disinfectants you use, particularly with respect to material and surface compatibility as well as the concentration information.
- ASSKEA GmbH recommends 'Sekusept® aktiv' for immersion disinfection of the accessories and 'Incidin® Plus' or 'Incidin® Liquid' for wipe disinfection of the device.

4.1.2 Cleaning and disinfecting the surface of the device

Clean the surfaces of the device regularly and disinfect them daily. All devices can be wiped with a damp, lint-free cloth. Follow the previous section 4.1.1 for wipe disinfection.

Repeated cleaning and disinfection procedures may result in minor discolorations of the plastic housing components. However these do not affect the function of the device.

If liquids or solids come into direct contact with the interior of the device, the device must be inspected by ASSKEA GmbH or by an authorized service partner of ASSKEA GmbH.

4.1.3 Cleaning and disinfecting the reusable secretion canister

Please note the respectively applicable hygiene instructions. Unless otherwise directed, please follow these steps:

1. Empty the reusable secretion canister and dispose of the aspirate in a proper manner.
2. Remove the ball before performing the following cleaning instructions (fig. 29(C)).
3. Rinse all components of the reusable secretion canister under running water.
4. Immerse all components of the reusable secretion canister in disinfectant solution for immersion disinfection (see section 4.1.1) recommended by ASSKEA GmbH, observing the specific concentration quantity. The reusable secretion canister and lid can also be autoclaved at 121°C for 20 minutes.


Fig. 29

- A** Reusable secretion canister lid
- B** Overflow protection
- C** Ball (for overflow protection)
- D** Reusable secretion canister

5. Then rinse the components thoroughly under running water and let them dry. When using an appropriate cleaner made by the manufacturer DR. WEIGERT (e.g. 'neodisher AN'), the components may also be cleaned in a special dishwasher, observing the manufacturer's instructions.
6. Assemble the dry components according to fig. 29.

ASSKEA GmbH recommends replacing the reusable secretion canister every 4 weeks, if it is used and disinfected frequently.

4.1.4 Cleaning and disinfecting the connecting tube

Please note the respectively applicable hygiene instructions. Unless otherwise directed, please follow these steps:

1. Applies only for ASSKEA reusable secretion canister system: Remove the external bacterial filter and dispose it of in a proper manner (see section 6.3) prior to cleaning the connecting tube.
2. Rinse the connecting tube with clear water and place it in disinfectant solution solution for immersion disinfection recommended by ASSKEA GmbH (see section 4.1.1), observing the specific concentration quantity and the residence times.
3. Then rinse the connecting tube thoroughly with clear water and let it dry.

ASSKEA GmbH recommends replacing the connecting tube every 4 weeks at the latest.

4.1.5 Cleaning and disinfecting the external canister 'Bag'

Please note the respectively applicable hygiene instructions. Unless otherwise directed, please follow these steps:

1. Rinse the external canister 'Bag' under running water.
2. Immerse the external canister 'Bag' in disinfectant solution for immersion disinfection (see section 4.1.1) recommended by ASSKEA GmbH, observing the specific concentration quantity.
3. Then rinse the external canister 'Bag' thoroughly under running water and let it dry.

You can also autoclave the external canister 'Bag' at 121°C for 20 minutes.

ASSKEA GmbH recommends replacing the external canister 'Bag' every 4 weeks at the latest, if it is used and disinfected frequently.

4.1.6 Cleaning/disposal of the suction tube

Please note the respectively applicable hygiene instructions. Unless otherwise directed, please follow these steps:

1. Dispose of the suction tube in a proper manner prior to each new patient!
2. If used for one patient only, rinse the suction tube with clear water (max. 0.5l) using the rinsing bottle after each aspiration procedure and place it in disinfectant solution for immersion disinfection recommended by ASSKEA GmbH (see section 4.1.1) daily, observing the specific concentration quantity and the residence times.
3. Then rinse the suction tube thoroughly with clear water and let it dry.

ASSKEA GmbH recommends replacing the suction tube for each new patient and no less than once per week at the latest, even if used for one patient only, since the material can become discolored and brittle due to frequent use.

4.1.7 Cleaning and disinfecting the rinsing bottle

Please note the respectively applicable hygiene instructions. Unless otherwise directed, please follow these steps:

1. Remove the lid of the rinsing bottle and rinse both parts under running water.
2. Immerse the rinsing bottle and the lid in disinfectant solution for immersion disinfection (see section 4.1.1) recommended

by ASSKEA GmbH, observing the specific concentration quantity.

3. Then rinse the components thoroughly under running water and let them dry.
4. Place the lid back on the rinsing bottle.

ASSKEA GmbH recommends replacing the rinsing bottle every 4 weeks, if it is used and disinfected frequently.

4.2 Cleaning and care in homecare use

4.2.1 General information

Health risks associated with handling infectious agents or pathogens

Infectious agents and pathogens in the aspirate cause health risks.

- Never operate the device without an external bacterial filter. (exception: ASSKEA disposable secretion canister system)
- Replace the external bacterial filter prior to each new patient or at least every two weeks when used exclusively for one patient, if it is used frequently! (exception: ASSKEA disposable secretion canister system)
- Reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!

- Replace the disposable liner 'One-Way' and suction tube, if it is full, prior to each new patient or no less than once per week.
- Prior to passing the device on to another patient or to other persons, professional hygienic processing by ASSKEA GmbH or an authorized service partner of ASSKEA GmbH is strictly required! For that purpose follow the information in section 4.3.
- After each aspiration procedure all components that have come into contact with the aspirate must be cleaned, disinfected or disposed of.
- Disposal of aspirate and contaminated components must be performed in a proper manner.

Health risks or damage to the device when handling disinfectants

If you have the opportunity to use disinfectants, please observe the information in section 4.1. Otherwise please follow the instructions in sections 4.2.2 through 4.2.7!

Potential bodily injury due to electric shock

- Switch off the device prior to cleaning/disinfection.
- Disconnect the power supply unit

or the mains cable by unplugging it from the mains supply. Then disconnect the power supply unit or the mains cable from the socket for power supply of the **ASSKEA secretion aspirator M- and S-series**.

4.2.2 Cleaning the surface of the device

Clean the surfaces of the device regularly. All devices can be wiped with a damp, lint-free cloth.

Repeated cleaning procedures may result in minor discolorations of the plastic housing components. However these do not affect the functioning of the device.

If liquids or solids come into direct contact with the interior of the device, the device must be inspected by ASSKEA GmbH or by an authorized service partner of ASSKEA GmbH.


4.2.3 Hygienic cleaning of the reusable secretion canister

1. Empty the reusable secretion canister and dispose of the aspirate in a proper manner.
2. Remove the ball before performing the following cleaning instructions (fig. 30(C)).
3. Rinse all components of the reusable se-

cretion canister under running water and let them dry.

4. Assemble the dry components according to fig. 30.

It is also possible to place all components of the

- A Reusable secretion canister lid
- B Overflow protection
- C Ball (for overflow protection)
- D Reusable secretion canister

Fig. 30

reusable secretion canister into a water bath at 65 °C for 15 minutes for hygienic cleaning. ASSKEA GmbH recommends replacing the

reusable secretion canister after 30 hygienic cleanings or in the case of visible damage or functional restrictions.

4.2.4 Hygienic cleaning of the connecting tube

1. Applies only to the ASSKEA reusable secretion canister system: Remove the external bacterial filter and dispose it of in a proper manner (see section 6.3) prior to cleaning the connecting tube.
2. Rinse the connecting tube with clear water and place it in a water bath at 65 °C for 15 minutes once a day.
3. Then rinse the connecting tube thoroughly with clear water and let it dry.

ASSKEA GmbH recommends replacing the connecting tube no less than every 4 weeks.

4.2.5 Hygienic cleaning of the external canister 'Bag'

Rinse the external canister 'Bag' under running water and let it dry.

It is also possible to place the external canister "Bag" into a water bath at 65 °C for 15 minutes for hygienic cleaning.

ASSKEA GmbH recommends replacing the external canister 'Bag' after 30 hygienic cleanings or in the case of visible damage or functional restrictions.

4.2.6 Hygienic cleaning of the suction tube

1. Dispose of the suction tube in a proper manner prior to each new patient!
2. If used for one patient only, rinse the suction tube with clear water (max. 0.5 l) using the rinsing bottle after each aspiration procedure and place it into a water bath at 65 °C for 15 minutes once a day.
3. Then rinse the suction tube thoroughly with clear water and let it dry.

ASSKEA GmbH recommends replacing the suction tube for each new patient and no less than every 4 weeks, even if used for one patient only, since the material can become discolored and brittle due to frequent use.

4.2.7 Hygienic cleaning of the rinsing bottle

1. Remove the lid of the rinsing bottle and rinse both parts under running water.
2. Let the components dry.
3. Place the lid back on the rinsing bottle.

It is also possible to place the rinsing bottle and the lid into a water bath at 65 °C for 15 minutes for hygienic cleaning.

ASSKEA GmbH recommends replacing the rinsing bottle after 30 hygienic cleanings or in the case of visible damage or functional restrictions.

4.3 Reuse of the device

The **ASSKEA secretion aspirators M- and S-series** are suitable for reuse. However, prior to passing the device on to other patients or persons, it must be processed professionally. Please send the **ASSKEA secretion aspirator M- and S-series** to ASSKEA GmbH or to qualified personnel authorized by ASSKEA GmbH. Follow the instructions in section 6.1!

4.4 Maintenance and service

Risk of performance impairments due to aging effects

If you notice impaired performance of the device or other problems, follow the instructions on troubleshooting in section 5 and contact the service, if necessary.

Risk associated with maintenance and service during operation

Maintenance and service must not be performed while the device is in use! The respective device must be switched off before maintenance and service work are performed.

Risk associated with performance of maintenance and service by unauthorized persons

Maintenance and service may be performed only by specialized personnel authorized by ASSKEA GmbH or by ASSKEA GmbH itself.

Risk associated with performance of maintenance and service without sufficient documents

Maintenance and service must be performed taking into account the relevant service instructions for the respective **ASSKEA secretion aspirator M- and S-series**. ASSKEA GmbH will also provide additional documents on request, if these can assist the service staff.

WARNING: Modification of the ME device is not permitted!

The **ASSKEA secretion aspirators M- and S-series** are maintenance-free, if used according to the instruction for use with the exception of components with limited shelf-life.

Perform a visual and functional check prior to each use. Make sure to include the accessories for the corresponding device.

Opening the device and repairs must be performed only by ASSKEA GmbH or by specialized personnel authorized by ASSKEA GmbH in compliance with the service documentation specified by the manufacturer as well as with technical and hygienic precautionary measures.

The device may be sent back for repair to ASSKEA GmbH directly or via the specialty dealer from whom you purchased the device.

Clean and disinfect all accessories prior to returning the device. The device itself must be treated with a surface disinfectant (see section 4.1.2 or 4.2.2). Please note the guidelines relating to decontamination before shipping (see section 6.1).

Please notify ASSKEA GmbH prior to returning the device to prevent time lags and label the shipping carton for the device with 'Used Medical Device'.

The corresponding product return form with the 'Used Medical Device' label can be found on the ASSKEA GmbH website at www.asskea.de, in the Service section under Product Return.

ASSKEA GmbH does not guarantee proper functioning of the **ASSKEA secretion aspirators M- and S-series** nor is ASSKEA GmbH liable for property damage or personal injury, if

- original ASSKEA accessories or spare parts are not used,
- the user instructions in this instruction for use are not followed,
- installation, settings, modifications, enhancements and/or repairs are not performed by ASSKEA GmbH or by personnel authorized by ASSKEA GmbH,
- the safety seal has been removed or is damaged.

4.5 Inspection of the ASSKEA secretion aspirators M- and S-series

ASSKEA GmbH offers its partners and customers fast and proper processing as well as required testing services.

5 Troubleshooting





5.1 ASSKEA secretion aspirators S-series

Error	Possible causes	Remedy
Device does not start.	Mains plug not plugged in	Check the electrical plug connection.
	No mains voltage	Please contact Service! (fuse defective)
	Device is too warm; temperature protection of the motor has been triggered.	Allow the device to cool down; temperature protection automatically resets after cooling.
insufficient flow rate	External bacterial filter polluted, moistened	Replace the external bacterial filter.
	Integrated bacterial filter polluted, moistened	Please contact Service!
	Leak in the tubing system	Check that secretion canister lid and tubing connections are secure and no leakages are apparent.
	Vacuum regulator (additional air valve) not closed	Close the vacuum regulator (additional air valve) (fig. 1(E)).
no flow rate	Manual overflow protection of the reusable secretion canister system is closed.	Check the fluid level in the secretion canister; empty the canister.
	Hydrophilic overflow protection of the disposable secretion canister system is blocked.	Replace the disposable liner 'OneWay'.
	External bacterial filter blocked	Replace the external bacterial filter.
	Integrated bacterial filter blocked	Please contact Service!
	Suction tube blocked	Clean/replace the suction tube.

Contact ASSKEA GmbH or your service partner if a malfunction cannot be corrected by the described measures.

5.2 ASSKEA secretion aspirators M-series

Error	LED display	Possible causes	Remedy
Device does not start.	 DFS  ON LEDs do not illuminate	Battery is discharged.	Connect the power supply unit for charging or operating.
	 LED illuminates red	Invalid input voltage	Connect the correct power supply unit; check the power supply.
Device turns itself off (applies to M20).	 LED blinks red, beeper on after 10s	Temperature of the battery too high	Switch off the device, if applicable, and allow it to cool down. Connect the power supply unit and switch on the device again.
	 ON LED switches off		
Motor turns itself off (applies to M20D, M30D).	 LED blinks red, beeper on after 10s, then the device turns itself off after a further 20s	Temperature of the battery too high	Switch off the device, if applicable, and allow it to cool down. Connect the power supply unit and switch on the device again.
	 ON LED switches off		
Device only runs with maximum vacuum (applies to M20D, M30D).	 LED illuminates red	Battery is discharged.	Switch off the device for charging after the aspiration procedure.
	 LED illuminates green when power supply unit is connected		
Device does not charge.	 LED blinks red	Temperature of the battery too high for charging	Allow the device to cool down. Automatic resumption of charging when battery has cooled down.
	 LED illuminates green		
	 LED illuminates red	Battery is defective.	Please contact Service!
	 LED illuminates green		
 LED illuminates red	Invalid input voltage	Connect the correct power supply unit, check the power supply.	

Error	LED display	Possible causes	Remedy
Insufficient flow rate	 LED illuminates red	Battery is weak.	Connect the power supply unit for charging or operating.
	 LED does not illuminate		
	 LED illuminates green/yellow	External bacterial filter blocked	Replace the external bacterial filter.
	 LED illuminates green	Integrated bacterial filter blocked	Please contact Service!
No flow rate	 LED illuminates green/yellow	Leak in the tubing system or in the canister system	Check that secretion canister lid and tubing connections are tight and no leakages are apparent.
	 LED illuminates green/yellow	Float ball of the overflow protection blocks the aspiration inlet.	Check the fluid level in the secretion canister, empty the canister.
	 LED illuminates green	Hydrophilic overflow protection of the disposable canister system is blocked.	Replace the disposable liner 'OneWay'.
	 LED illuminates green	External bacterial filter blocked	Replace the external bacterial filter.
	 LED illuminates red (only for M20D and M30D)	Integrated bacterial filter blocked	Please contact Service!
	 DFS	Suction tube blocked	Clean/replace the suction tube.

Contact ASSKEA GmbH or your service partner if a malfunction cannot be corrected by the described measures.

6 Transport, storage and disposal

6.1 Decontamination prior to shipment

Prior to transferring the **ASSKEA secretion aspirators M- and S-series** to another patient, the devices must be properly processed by ASSKEA GmbH or by specialized personnel authorized by ASSKEA GmbH to protect subsequent users and other persons. Processing is mandatory in accordance with the German Medical Device Operator Ordinance (MPBetreibV), German Medical Devices Act (MPG) and the manufacturer's instructions.

ASSKEA GmbH offers its partners and customers fast and proper processing as well as required testing services.

The **ASSKEA secretion aspirators M- and S-series** must be cleaned and disinfected prior to shipment to ASSKEA GmbH. Please follow the instructions in section 4.1.2 or 4.2.2! Affix the 'Used Medical Device' label provided to the shipping carton! Please give ASSKEA GmbH advance notice of your product return. The product return form can be found on our website at www.asskea.de in the Service section under Product Return.

6.2 Storage

Store the **ASSKEA secretion aspirators M- and S-series** as indicated in the technical data (see section 7)! Store the device in the shipping carton till next use.

The battery of the **ASSKEA secretion aspirators M-series** must be charged prior to storage of the device. This ensures that the device is operational at all times.

Fully recharge the battery if the **ASSKEA secretion aspirator M-series** is not used for a prolonged period of time (approx. 10 months)!

6.3 Disposal

- Disposal of the device and of the accessories must be carried out in a proper manner.
- Decontaminate the device and the accessories prior to disposal.
- According to EU Directives 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) and 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II), the device must not be disposed of in domestic waste.
- The devices are registered in the national register for waste electric equipment (EAR foundation) as small electrical devices and can be submitted to a collection facility nearby for disposal.
- The device and accessories may also be disposed of via ASSKEA GmbH or the service partner.
- Outside of the EU: Follow the disposal requirements of your country!

7 Technical data

7.1 ASSKEA secretion aspirators S-series

Flow rate* (measuring point at suction tube nozzle)	ASSKEA S20: 18 l/min ± 2 l/min (low flow) ASSKEA S20K: 18 l/min ± 3 l/min (low flow) ASSKEA S30: 28 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg) (high vacuum) conversion factor: 1 kPa ~ 7.5 mmHg; 1 kPa ~ 10 mbar
Canister	Reusable secretion canister system (1,000 ml); Disposable secretion canister system (1,000 ml)
Suction tube (item number: 100713-2)	PVC, inner diameter 6 mm, wall thickness 2 mm, length 150 cm, vacuum stable up to -0.9 bar
Mains voltage	230 V AC
Mains frequency	50 Hz
Power consumption	ASSKEA S20: 80 W ASSKEA S20K: 156 W ASSKEA S30: 80 W
Current consumption at 230 V	ASSKEA S20: 0.35 A ASSKEA S20K 0.72 A ASSKEA S30: 0.35 A
Dimensions of the basic device (H x W x D)	290 mm x 259 mm x 130 mm
Weight (basic device)	ASSKEA S20: 3.4 kg ASSKEA S20K: 3.5 kg ASSKEA S30: 4.1 kg
Operating time	20 minutes ON, 40 minutes OFF
Service life	5 years

Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pursuant to 93/42/EEC, IX	Ila
Protection class pursuant to IEC 60601-1	II
IP degree of protection	IP22
CE marking	CE1434
Sound emission	without vacuum: ASSKEA S20: approx. 55 dB(A) ASSKEA S20K: approx. 66 dB(A) ASSKEA S30: approx. 57 dB(A)
Ambient conditions	Transport/storage: -20 °C to +60 °C Operation: 0 °C to +35 °C Relative humidity: 5% to 80%, non-condensing Air pressure: 800 hPa to 1060 hPa (at maximum 2000 m altitude above sea level)
Mains cable	H05VV-F 2x1,0 mm ² , length: 2.5 m
Item number (REF)	ASSKEA S20: 100009-3 ASSKEA S20K: 100011-3 ASSKEA S30: 100010-3

* The information provided may differ depending on the altitude above sea level, the prevailing air pressure and the air temperature.

7.2 ASSKEA secretion aspirators M-series

Flow rate* (measuring point at suction tube nozzle)	ASSKEA M20: 18 l/min ± 2 l/min (low flow) ASSKEA M20D: 18 l/min ± 2 l/min (low flow) ASSKEA M30D: 26 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg) (high vacuum) conversion factor: 1 kPa ~ 7.5 mmHg; 1kPa ~ 10 mbar
Canister	Reusable secretion canister system (1,000 ml); Disposable secretion canister system (1,000 ml)
Suction tube (item number: 100713-2)	PVC, inner diameter 6 mm, wall thickness 2 mm, length 150 cm, vacuum stable up to -0.9 bar
Nominal voltage of the power supply unit	In: AC 100 – 240 V~/ 50-60Hz / 1.5 A Out: DC 12V / 5A
Maximum load current	5.0 A
Permissible input voltage	12 V
Power consumption at 12 V	ASSKEA M20: 47 W ASSKEA M20D: 57 W ASSKEA M30D: 52 W
Power input at 12 V	ASSKEA M20: 3.5 A ASSKEA M20D: 3.5 A ASSKEA M30D: 3.8 A
Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pursuant to 93/42/EEC, IX	IIa
Protection class pursuant to IEC 60601-1	II
IP degree of protection	IP22
CE marking	CE1434
UL marking	E355754
Sound emission	without vacuum: approx. 62 dB(A)

Ambient conditions	Transport/storage: -20 °C to +60 °C Operation: 0 °C to +35 °C Relative humidity: 5 % to 80 %, non-condensing Air pressure: 800 hPa to 1060 hPa (at maximum 2000 m altitude above sea level)
Battery, rechargeable	10.8 V; 5800 mAh (lithium ion battery)
Charging time, if battery is empty	2.5 h – 3 h
Charging time, if battery is approx. 50 % full	1.5 h – 2 h
Power supply unit	GTM91099-6015-3.0-T2
Mains cable of the power supply unit	H03VVH2-Fx0,75 2.5 A 250 V~; length: 1.8 m
Dimensions of the basic device (H x W x D)	290 mm x 259 mm x 130 mm
Weight (basic device)	ASSKEA M20: 3.1 kg ASSKEA M20D: 3.1 kg ASSKEA M30D: 3.7 kg
Operating time	ASSKEA M20: approx. 45 minutes ON, 10 minutes OFF ASSKEA M20D: approx. 45 minutes ON, 10 minutes OFF ASSKEA M30D: approx. 40 minutes ON, 10 minutes OFF
Service life	5 years
Runtime in battery operation	ASSKEA M20: approx. 45 minutes ASSKEA M20D: approx. 45 minutes ASSKEA M30D: approx. 40 minutes
Item number (REF)	ASSKEA M20: 100003-3 ASSKEA M20D: 100005-3 ASSKEA M30D: 100006-3

* The information provided may differ depending on the altitude above sea level, the prevailing air pressure and the air temperature.

8 EMC information

WARNING: Portable and mobile RF communication equipment (incl. peripheral devices such as antenna cables and external antennas) may influence medical electrical devices and therefore should not be used within a range of 30 cm of any part of the **ASSKEA secretion aspirators M- and S-series** including the cables. Otherwise the performance of the device may be impaired.

WARNING: Use of accessories and spare parts, transformers and cables for the **ASSKEA secretion aspirators M- and S-series** not indicated or provided by ASSKEA GmbH may increase the electromagnetic emissions or reduce the electromagnetic immunity of the devices, resulting in impermissible operation.

No warranty is provided for damage caused by using accessories and spare parts, transformers and cables not recommended or by improper use. Use only original ASSKEA accessories and spare parts!

WARNING: Use of the indicated or provided accessories and spare parts, transformers and cables for other devices than the **ASSKEA secretion aspirators M- and S-series** may increase the electromagnetic emissions or reduce the electromagnetic immunity. No warranty is provided for damage caused by using the indicated or provided accessories and spare parts, transformers and cables with other devices or by improper use. Use only the accessories and spare parts, transformers and cables with the **ASSKEA secretion aspirators M- and S-series!**

WARNING: Use of the **ASSKEA secretion aspirators M- and S-series** directly adjacent to or stacked with other devices should be avoided, since this could lead to impermissible operation. If use of the **ASSKEA secretion aspirators M- and S-series** adjacent to or stacked with other devices is required, the **ASSKEA secretion aspirators M- and S-series** and the other devices should be monitored in order to verify proper operation in this arrangement!

WARNING: The **ASSKEA secretion aspirators M- and S-series** may influence other devices, examinations and treatments electromagnetically. For this reason, special attention should always be paid to other devices and examinations or treatments performed in parallel so that any influence is detected as soon as possible.

The **ASSKEA secretion aspirators M- and S-series** meet the requirements of IEC 60601-1-2/EN 60601-1-2 'Electromagnetic Compatibility – Medical Electrical Equipment' without deviations or restrictions. Therefore electromagnetic interference is reduced to a minimum.

Follow the indicated instructions and guidelines to sustain the basic safety and the essential functions of the **ASSKEA secretion aspirators M- and S-series** over the entire expected service life.

8.1 Electromagnetic environment, in which the ASSKEA secretion aspirators M- and S-series may be operated

The **ASSKEA secretion aspirators M- and S-series** are intended for operation in the electromagnetic environment specified below, in which RF disturbances are controlled. The customer or the user of the **ASSKEA secretion aspirators M- and S-series** must ensure that they are operated in such an environment.

The environments for intended operation include professional healthcare institutions and the home care environment. Special environments, such as close to RF surgery or MRI or environments, in which the intensity of EMC disturbances is high, are excluded.

Emission limits

Conducted and radiated RF emissions	CISPR 11
-------------------------------------	----------

Housing		
Phenomenon	Test method	Immunity test levels
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 15 kV air
Radiated RF disturbances	IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz
Radiated RF disturbances	pursuant to frequencies and immunity test levels of EN 60601-1-2, table 9	see EN 60601-1-2, Table 9
Magnetic field at rated power frequency (50 Hz)	IEC 61000-4-8	30 A/m

Power supply AC		
Phenomenon	Test method	Immunity test levels
Fast transient electrical disturbances/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency
Impulse voltages/surges	IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground
Conducted RF disturbances	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Voltage dips	IEC 61000-4-11	0 % U_T for ½ period At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0 % U_T for 1 period 70 % U_T for 25 periods one-phase at 0°
Voltage interruptions	IEC 61000-4-11	0 % U_T for 250 cycles

Note: U_T is the AC mains voltage prior to the application of the test levels.

8.2 Dealing with electromagnetic interactions

Although electromagnetic interference of the **ASSKEA secretion aspirators M- and S-series** has been reduced to a minimum, electromagnetic disturbances between the **ASSKEA secretion aspirators M- and S-series** and other devices cannot be excluded. Therefore you should comply with the specified requirements and instructions regarding the permissible electromagnetic environment in any case and monitor the **ASSKEA secretion aspirators M- and S-series** in order to ensure proper operation and to prevent adverse events for patients and users. Select another location for the **ASSKEA secretion aspirators M- and S-series** if the permissible electromagnetic environment cannot be ensured or if you have observed functioning of the **ASSKEA secretion aspirators M- and S-series** or of other devices in the vicinity that is not intended.

Since electromagnetic propagation is affected by absorption and reflection from structures, objects and people and since the field strength from fixed transmitters cannot be predicted in advance with accuracy, a site survey of the electromagnetic phenomena at the location at which the **ASSKEA secretion aspirators M- and S-series** shall be operated should be considered to assess the existing electromagnetic environment at this location. If unusual performance is detected, additional measures may be required, for example a modified orientation or another location for the **ASSKEA secretion aspirators M- and S-series**. If the essential functions of the **ASSKEA secretion aspirators M- and S-series** are impaired by electromagnetic disturbances, a reduced or non-existent flow rate and vacuum may be expected.

8.3 List of all cables and transformers replaceable by the user

a) ASSKEA S-series

Name	Specification	Maximum length
Mains cable	Type: H05VV-F 2x1,0 mm ² Technical data: 230 V~, 50 Hz	2.5 m

b) ASSKEA M-series

Name	Specification	Maximum length
Power supply unit incl. connecting cable	Type: GTM91099-6015-3.0-T2 Technical data: 100-240 V~, 50-60 Hz, 1.5 A (in) 12 V DC, 5 A (out)	1.20 m
Mains cable	Type: H03VVH2-F Technical data: 250 V~, 2.5 A	1.80 m

9 Order information

Item number	Description	PU
100000	External canister 'Bag'	1
100002	Disposable liner 'OneWay' (1,000 ml)	60
100267	Holder for external canister 'Bag'	1
100013	Double filter system replacement set (DFS®)	1
100280	Connecting tube for disposable secretion canister system	1
100288	Universal bed holder	1
100346	Variable holder for pipe and rail systems	1
100385	Disposable secretion canister system	1
100295	Bag for ASSKEA aspirators	1
100194	External bacterial filter	5
100278	Ball for overflow protection	3
100279	Reusable secretion canister (1,000 ml)	1
100205	S6 lid for reusable secretion canister	1
100199	S6 reusable secretion canister complete	1
100282	Rinsing bottle (250 ml)	1
100283	Tube mount (holding ring with tube clamp)	1
100759	Power supply unit GTM91099-6015-3.0-T2 for ASSKEA M-series incl. mains cable	1
100216-1	Mains cable for ASSKEA S-series	1
100753	Car connecting cable for ASSKEA M-series	1
100212-2	S6 tubing set (Ø 6 mm)	1
100713-2	Suction tube with funnels and fingertip (unsterile)	10
100294	S6 tube connector	5
100222	S6 fingertip	5
100202-2	S6 reusable secretion canister system (incl. tubing set)	1
100732-1	Stand with 5 rollers	1
100736	Device holder for roller stand	1
100704	Wall holder stainless steel	1

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Dok.-Nr.: GA01-02 m, Revision J, 12.2019

CE 1434

MEDICAL — GENERAL MEDICAL EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012);
CAN/CSA C22.2 No. 60601-1:14